wallace[®]

Malleable Stylets

en

INSTRUCTIONS FOR USE

fr

MODE D'EMPLOI

de

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it

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es

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pt

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nl

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sv

BRUKSANVISNING

da

BRUGSANVISNING

no

BRUKSANVISNING

fi

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el

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Wallace® Malleable Stylets

INSTRUCTIONS FOR USE

Wallace® malleable stylets are closed ended single-use devices provided sterile for use with Wallace® embryo transfer catheters, should passage through the cervix be impeded.

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

WARNING:

- **THESE PROCEDURES SHOULD ONLY BE PERFORMED BY PERSONS HAVING ADEQUATE TRAINING AND FAMILIARITY WITH THESE TECHNIQUES. CONSULT MEDICAL LITERATURE REGARDING TECHNIQUES, COMPLICATIONS AND HAZARDS PRIOR TO PERFORMANCE OF THESE PROCEDURES.**
- **TO BE USED BY, OR UNDER THE DIRECTION OF, QUALIFIED PERSONS IN LINE WITH LOCAL GUIDELINES GOVERNING IN VITRO FERTILISATION, IF APPLICABLE.**
- **DISCARD IF PRODUCT OR PACKAGING IS DAMAGED.**

Description

These instructions are intended for use with sterile Wallace® malleable stylets product codes

Length	For use with	Code
18cm	1816, 1816S, TT1816, CE118, CE418	1816ST
23cm	1816N, TT1816N, 2316, CE123, CE423	1816NST

The devices consist of a malleable inner stylet and a detachable outer sheath that is attached to the inner stylet by a Luer compatible adapter. With the hubs locked together the inner malleable stylet forms a smooth radius at the tip of the outer sheath. The outer sheath has a series of 1cm graduations at the distal end.

Indications for use

To provide rigidity when using a Wallace® Embryo Replacement Catheter and passage through the cervix is impeded or difficult.

Contraindications

The catheter should not be used:

- In the presence of chronic cervical infection
- In the presence of or after recent pelvic inflammatory disease

Instructions for Use

Place patient in the lithotomy, dorsal, or left lateral position. The cervix should be exposed with a speculum and gently swabbed with cotton wool moistened with normal saline or medium.

With the hubs locked in position, advance the malleable stylet through the cervical canal to the internal os.

WARNING

The catheter should never be forced against digitally felt resistance while inside the uterine cavity, as forcing the catheter may result in damage to the endometrial tissue and bleeding.

The outer sheath should not be advanced further than the internal os, and should certainly never enter the uterine cavity, as this may result in damage to the endometrial tissue and bleeding.

The inner and outer hubs must be locked together. Removing the stylet from the outer sheath, before placement, may result in damage to the endometrial tissue and bleeding.

If slight resistance is experienced, withdraw the catheter. Using thumb and forefinger curve the malleable stylet to complement the patient's anatomy and pass the catheter in this position to the internal os.

Gently uncouple and remove the inner stylet leaving the outer sheath in place. Complete the embryo replacement in accordance with Wallace® embryo replacement Instructions for use.

Dispose of the stylet in accordance with local medical hazardous waste practices.

The above are only suggested instructions for use. The final decision on the technique used is the responsibility of the clinician in charge.

Wallace, Wallace design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2014 Smiths Medical. All rights reserved.

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Wallace, designmerket Wallace og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet ® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.
Wallace, znak towarowy Wallace oraz znak towarowy Smiths Medical są znakami handlowymi Smiths Medical. Symbol® oznacza, że znak towarowy jest zarejestrowany w Urzędzie Patentowym Stanów Zjednoczonych (U.S. Patent and Trademark Office) i w urzędach patentowych innych krajów. Wszystkie inne wymienione nazwy i znaki są nazwami towarowymi, znakami towarowymi lub nazwami usług ich właścicieli.