***Nikon***

**Upright Microscope**

**ECLIPSE**

***Si***

**Instructions**

## Introduction

Thank you for purchasing a Nikon product.

This instruction manual is intended for users of the Nikon Upright Microscope ECLIPSE Si RS.

No part of this manual may be reproduced or transmitted in any form without prior written permission from Nikon.

- The contents of this manual are subject to change without notice.
  - The equipment described in this manual might differ from the actual product in its appearance.
  - Although every effort has been made to ensure the accuracy of this manual, errors or inconsistencies might remain.
  - If you notice any points that are unclear or incorrect, please contact your local Nikon representative.
  - Some of the equipment described in this manual may not be included in the set you have purchased.
  - If you intend to use any other equipment with this product, read the manual for that equipment too.
  - If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment might be impaired.
-

## Contents

Safety Precautions .....	3
Meaning of Symbols Used on the Product.....	3
WARNING and CAUTION Symbols .....	3
⚠ WARNING .....	4
⚠ CAUTION .....	6
Sicherheitsvorkehrungen.....	7
Consignes de sécurité .....	12
Precauciones de seguridad .....	17
Precauzioni per la sicurezza .....	22
Veiligheidsvoorschriften.....	27
Sikkerhedsforanstaltninger .....	32
Precauções de segurança .....	37
Turvatoimet.....	42
Προφυλάξεις ασφαλείας.....	46
Notes on Handling the Product .....	51
<b>Chapter 1 Nomenclature of Each Part .....</b>	<b>54</b>
<b>Chapter 2 Microscopy Techniques .....</b>	<b>56</b>
2.1 Bright-field Microscopy .....	57
2.2 Phase Contrast Microscopy .....	64
2.3 Simple Polarizing Microscopy .....	68
2.4 Diascopic Fluorescence Microscopy.....	72
2.5 Dark-field Microscopy.....	76
<b>Chapter 3 Detailed Explanation.....</b>	<b>79</b>
3.1 Turning the Power On and Off .....	79
3.1.1 Turn On/Off the Power to the Microscope .....	79
3.1.2 ECO Mode and SLEEP Mode .....	79
3.2 Operating Diascopic Illumination.....	80
3.2.1 Turning Diascopic Illumination On and Off and Adjusting the Brightness .....	80
3.2.2 Using Light Intensity Management .....	80
3.3 Operating the Stage Specimen Holder .....	81
3.3.1 Placing a Specimen Using the Specimen Holder .....	81
3.3.2 Operating the Stage.....	81
3.4 Operating the Nosepiece and Objectives .....	82
3.4.1 Switching Objectives.....	82
3.4.2 Oil-Immersion Observation .....	82
3.4.3 Performing Water-Immersion Microscopy .....	83
3.5 Focusing Operations .....	84
3.5.1 Focusing on a Specimen .....	84
3.5.2 Torque Adjustment of the Coarse Focus Knob .....	85
3.5.3 Replacing a Specimen Using the Stage Vertical Movement Stopper.....	86
3.6 Operating the Binoculars.....	87
3.6.1 Adjustment of Diopter .....	87
3.6.2 Adjusting Interpupillary Distance .....	88
3.7 Operating the Condenser.....	89
3.7.1 Focusing the Condenser.....	89
3.7.2 Centering the Condenser.....	89
3.8 Adjustment of Aperture Diaphragm.....	90
3.9 Adjustment of field diaphragm.....	91

3.10	Tips for Each Microscopy Method .....	92
3.10.1	Tips for Phase Contrast Microscopy .....	92
3.10.2	Tips for Simple Polarizing Microscopy .....	93
3.10.3	Tips for Diascopic Fluorescence Microscopy .....	94
3.10.4	Tips for Dark-field Microscopy .....	95
3.11	Combinations of Optical Parts.....	96
3.12	Microscope Terminology .....	97
<b>Chapter 4</b>	<b>Assembly .....</b>	<b>98</b>
4.1	Configuration Diagrams .....	98
4.2	Assembly.....	99
4.2.1	Basic Configuration (for Bright-field Microscopy) .....	99
4.2.2	Configuration for Phase Contrast Microscopy .....	103
4.2.3	Configuration for Simple Polarizing Microscopy .....	104
4.2.4	Configuration for Diascopic Fluorescence Microscopy.....	105
4.2.5	Configuration for Dark-Field Microscopy .....	108
4.3	Setting Up the Microscope .....	109
<b>Chapter 5</b>	<b>Troubleshooting .....</b>	<b>111</b>
5.1	Problems With the Microscope Image .....	111
5.1.1	Common Problems for Each Microscopy .....	111
5.1.2	Problems in Phase Contrast Microscopy.....	112
5.1.3	Problems in Simple Polarizing Microscopy.....	113
5.1.4	Problems in Diascopic Fluorescence Microscopy .....	113
5.1.5	Problems in Dark-field Microscopy .....	114
5.2	Problems in Operation .....	115
5.3	Electrical.....	116
<b>Chapter 6</b>	<b>Care and Maintenance .....</b>	<b>117</b>
6.1	Cleaning the Lenses .....	117
6.2	Cleaning the Microscope.....	117
6.3	Disinfecting the Microscope .....	117
6.4	Storing the Microscope .....	117
6.5	Periodical Inspections (Charged).....	117
<b>Chapter 7</b>	<b>Specifications .....</b>	<b>118</b>
7.1	Microscopy (Principles).....	118
7.2	Performance Properties .....	119
7.3	Physical Properties .....	120
7.4	AC adapter .....	120
7.5	Power Cord .....	121
7.6	Safety standards .....	121
7.7	External dimensions (unit: mm).....	122


## Safety Precautions

Although this product is designed and manufactured to be completely safe during use, incorrect usage or failure to follow the safety instructions provided may cause personal injury or property damage. To ensure correct usage, read this manual carefully before using the product. Do not discard this manual and keep it handy for easy reference.

### Meaning of Symbols Used on the Product

Symbols that appear on the product indicate the need for caution at all times during use. Always refer to the instruction manual and read the relevant instructions before manipulating any part to which the symbol has been affixed.

Symbol	Contents
	<p><b>Biohazard</b></p> <p>This symbol can be found on the stage of the microscope body (stage handle side), and cautions the following:</p> <ul style="list-style-type: none"> <li>• <b>WARNING:</b> The product may become biohazardous if a sample is spilled onto the product.</li> <li>• To avoid exposure to biohazard, do not touch contaminated parts with your bare hands.</li> <li>• Decontaminate the contaminated parts according to the standard procedures for your facility.</li> </ul>
	<p><b>WARNING (photobiological safety)</b></p> <p>This symbol can be found on the microscope body (stage handle side), and cautions the following:</p> <ul style="list-style-type: none"> <li>• Photobiological safety risk group classification of diascope illumination</li> <li>• Do not look into diascope illumination.</li> </ul> <p>For details, see  <b>WARNING</b>, “7. Photobiological safety” and “8. Do not look into the illumination section.”</p>

Microscope main body (stage handle side)

### WARNING and CAUTION Symbols

Safety instructions in this manual are marked with the following symbols to highlight their importance. For your safety, always follow the instructions marked with these symbols.

Symbol	Contents
 <b>WARNING</b>	Disregarding instructions marked with this symbol may lead to serious injury or death.
 <b>CAUTION</b>	Disregarding instructions marked with this symbol may lead to injury or property damage.

## WARNING

### 1. Do not disassemble.

Disassembling this product may result in electric shock or malfunction. Malfunctions and damage due to disassembly will not be warranted. Do not disassemble any part unless instructed to do so in this manual. If you experience problems with the product, contact your nearest Nikon representative.

### 2. Read the instructions thoroughly.

To ensure safety, thoroughly read this manual and the manuals for other equipment to be used with this product. In particular, be sure to follow the warnings and cautions at the beginning of the manuals.

### 3. Input voltage, AC adapter, power cord

The AC adapter for supplying power to this product can be used with 100 to 240 VAC at 50 to 60 Hz. Use only the AC adapter and power cord specified in "Chapter 7 Specifications." Use of another AC adapter or power cord may result in malfunction or fire.

- Avoid using the product in an environment where the supply voltage may fluctuate excessively.
- The AC adapter for this product is classified as Class I for electrical shock protection. Be sure to connect it to a protective earth terminal.
- If the AC adapter is covered or items are placed on the AC adapter, heat dissipation may be hindered, causing it to become abnormally hot.

### 4. Handling flammable solvents

The following flammable solvents are used with the product:

- Immersion oil (Nikon Immersion Oil for oil immersion objectives)
- Absolute alcohol (ethyl alcohol or methyl alcohol for cleaning optical parts)
- Petroleum benzine (for wiping off the immersion oil)
- Medical alcohol (for disinfecting the microscope)

Never hold a flame near these solvents. When using a solvent, ensure correct and safe use by thoroughly reading the instructions provided by the solvent manufacturer.

When using a solvent with this product, observe the following precautions:

- Keep solvents away from the product and its surroundings when plugging/unplugging the AC adapter or power cord.
- Be careful not to spill the solvents.

### 5. Hazardous samples

This microscope is mainly for use in microscopic observation of a sample such as cells and tissue affixed to a slide.

When handling a specimen, check to determine whether the specimen is hazardous. Handle hazardous specimens according to the standard procedure for your laboratory. If the specimen is of an infectious nature, wear rubber gloves to avoid infection, and be careful not to touch a specimen.

If such a sample comes in contact with the microscope, decontaminate the contaminated portion according to the standard procedure for your laboratory.

### 6. Notes on operating the sliders on the condenser part

A slider for microscopy has two set positions. It is attached to the condenser part and used for switching the observation methods or magnifications by changing the set position. The bright light might reach the eyepiece tube when you move the slider to change the set position.

Remove your eyes from the binoculars or turn off the diasopic illumination when you operate a slider.

Especially in diasopic fluorescence microscopy, the very bright illumination is used, so the intense flashes might reach the eyepiece tube when you operate the slider.

## WARNING

### 7. Photobiological safety

This product is designed and manufactured in accordance with the safety standard IEC 62471 “Photobiological safety of lamps and lamp systems.”

Illumination near the stage and light from the binocular part of the eyepiece tube or the teaching head, and light from the camera port of the trinocular eyepiece tube are classified into the following risk group. The distance (hazard distance) from the vicinity of the stage or from the binocular part of the eyepiece tube or the teaching head, or from the opening of the camera port where the risk group classification is equivalent to the exempt group that does not cause photobiological injury is as follows.

	Risk group classification	Hazard distance
<b>Retinal blue light hazard</b>	Risk group 3	36 m
<b>Infrared radiation hazard to eyes</b>	Risk group 3	127 m

Illumination in risk group 3 can be hazardous. Pay attention to brightness adjustment by the brightness control knob and avoid looking strong light for a long time.

- The LED light source is built in at the bottom part of the microscope. Illumination light from the field lens goes through the condenser and illuminates the specimen. Do not stare at the condenser lens, specimen or periphery of the specimen, and do not look into the field lens.
- Do not look into strong light from the binocular part of the eyepiece tube or the teaching head.
- Do not remove the eyepieces from the binocular part when illumination is turned on.
- Do not look into strong light from the camera port from directly above.
- When you do not attach the camera, attach the supplied cap to the camera port.

### 8. Do not look into the illumination section.

The following warning label indicating the highest risk (Risk Group 3) to photobiological safety is affixed on the handle side of the stage to remind the user of the following precautions. (For the locations of this label, see the figure in “Chapter 1 Nomenclature of Each Part.”)

#### **WARNING (Risk Group 3)**

IR and possibly hazardous light is emitted from the camera port when the trinocular eyepiece tube is used.

Never look into the camera port.

#### **CAUTION (Risk Group 2)**

Possibly hazardous light is emitted from the condenser lens when diasopic illumination is used.

Do not stare at emitted light. It may be harmful to the eyes.

#### **NOTICE (Risk Group 1)**

Possibly hazardous light is emitted from the eyepiece lenses when diasopic illumination is used.

Avoid long-time microscopic observation with strong illumination.

## CAUTION

### 1. Assembly of the microscope

- Assemble the microscope while the AC adapter is not connected to it.
- Take care to avoid pinching your fingers and hands.

### 2. Do not wet the product or allow ingress of foreign matter.



Do not allow the product to become wet, as doing so may result in malfunction, overheating, or electric shock.

If water or other liquids are accidentally spilled onto the product, immediately unplug the power cord from the AC adapter. Then, wipe off the liquid with a dry cloth.


Ingress of foreign matter may also result in malfunction. If liquids or foreign matter gets in the product, cease use of the product, and contact your nearest Nikon representative.

### 3. Moving the product

- The microscope has carrying handholds and handles.  
There are two handholds on both sides of the base and two handles on the arm on the rear side of the microscope.  
When carrying the microscope, hold the microscope firmly by these handholds and handles.
- The storage part on the rear of the microscope is not a handle. Do not hold this part as it stores the AC adapter and power cord.
- When moving the product, do not hold by the focus knobs, eyepiece tube, and stage, etc. The parts may become detached and also result in malfunctions.

Hold a handle and a handhold.

Hold the handles.

### 4. Disposal of the microscope

To avoid biohazard risks, dispose of the microscope as contaminated equipment, according to the standard procedures for your facility.

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## Notes on Handling the Product

### 1. Handle with care.

The product is a precision optical instrument. Handle the product with care and avoid physical shocks and vibrations.

In particular, the precision of objectives may be lost by even weak physical shocks.

### 2. Electromagnetic environment

Before using this product, Nikon recommends evaluating the electromagnetic environment of the installation site.

Do not use this product close to strong electromagnetic radiation sources (example: unshielded intentional RF sources.) They may interfere with the proper operation of this product.

This product emits low-level electromagnetic radiation. Do not install this product near precision electronic devices. Otherwise, the performance of such devices might be degraded. If TV or radio reception is affected, move the TV or radio farther away from this product.

This product complies with requirements of EMC directive (2014/30/EU) and IVD directive (98/79/EC).

- Achieve electromagnetic compatibility to ensure satisfactory operation of the product.  
This product is tested and in conformity with FCC regulations (FCC Part15 Subpart B: Radio Frequency Device (commercial and industrial areas)).
- Use this product only in the commercial or industrial area. Using this product in the residential area may affect other devices.

### 3. Installation location and storage location

The product is a precision optical instrument. Use or storage of the product under inappropriate conditions may result in malfunctions or loss of precision. Details on the installation and storage environment are given in Chapter 7, "Specifications." Additionally, the following conditions must be considered for the installation location and the storage location.


- Avoid placing the product in a hot or humid location. Use or storage of the product in a hot or humid location may result in molding of or condensation on the lenses, loss of precision, and malfunctions.
- Install the product in a place with little dust and dirt.  
When storing the product, place a cover over the product to protect it from dust.
- Install the product in a place with little vibration.
- Install and store the product on a level and sturdy table or stage that can bear the weight of the product.  
Install the product in a location with minimal exposure to hazards in the event of earthquakes and other potential disasters. If necessary, secure the product to the working desk or other heavy and stable items with a strong wire or other means, so as to prevent it from falling.
- Install the product so that the power cord can be unplugged immediately from the AC adapter in case of an emergency.
- Avoid placing the product in direct sunlight or immediately under room lights.  
The image quality is degraded in a bright environment due to the extraneous light entering the objective. A room light immediately above the microscope may also enter the objective as extraneous light. In this case, it is recommended that you turn off the room light above the microscope.
- Install the product at least 10 cm away from the surrounding walls.
- Do not place items on the product.

### 4. Handling optical parts

Handle optical components such as lenses and filters with care, so as not to damage them. Scratches and dirt (i.e., fingerprints) on the optical components will degrade the microscope image. If components require cleaning, clean them by following the procedure in "Chapter 6 Care and Maintenance."

## 5. Fastening plate

The microscope is held tightly by the red fastening plate during shipment. Be sure to remove the plate before use. The microscope cannot be used without removing the fastening plate. For details, see “Chapter 4 Assembly.”

---

This page is intentionally blank.



---

# 1


## Nomenclature of Each Part

This chapter describes the names of parts and the functions of controls of the microscope system.


### Nomenclature of Each Part


Lower side of the binoculars

Eyepiece tube mount


Rear of the microscope

AC adapter and power cord

## Functions of Controls

- |  |  |
|--|--|
| <p><b>(1) Diopter ring</b><br/>Compensates for the difference between the right and left eyesight of the observer.</p> <p><b>(2) Nosepiece</b><br/>Switches the objectives.</p> <p><b>(3) Specimen holder</b><br/>Holds a specimen (glass slide).</p> <p><b>(4) Stage movement knob</b><br/>Moves the specimen holder back and forth (upper knob) and right and left (lower knob).</p> <p><b>(5) Stage vertical movement stopper</b><br/>Sets the upper limit of the stage movement. Used for the fail-safe device for specimen and objective damage protection.</p> <p><b>(6) Aperture diaphragm lever</b><br/>Adjusts the aperture diaphragm from <math>\phi 1.5</math> mm to <math>\phi 18</math> mm.</p> <p><b>(7) Condenser vertical movement knob</b><br/>Vertically moves the condenser (focusing).</p> <p><b>(8) Condenser centering knob</b><br/>Horizontally moves the condenser (centering).</p> <p><b>(9) Coarse focus knob torque adjustment ring</b><br/>Adjusts the tension of the coarse focus knob.</p> | <p><b>(10) Coarse focus knob</b><br/>Used for coarse stage movement (coarse focusing).</p> <p><b>(10) Fine focus knob</b><br/>Used for fine stage movement (fine focusing).</p> <p><b>(12) Field diaphragm lever</b><br/>Adjusts the field diaphragm from <math>\phi 1.5</math> mm to <math>\phi 28</math> mm.</p> <p><b>(13) Status display</b><br/>Displays the nosepiece address, objective name, ECO mode state, LIM on/off state, illumination on/off state and brightness control status bar.</p> <p><b>(14) Illumination on/off button</b><br/>Short-press: Turns illumination on or off.<br/>Long-press: Activates the SLEEP mode.</p> <p><b>(15) Brightness control knob</b><br/><b>Controls diasopic illumination.</b><br/><b>Turn:</b> Adjust the illumination intensity. (clockwise: brighter, counterclockwise: darker)<br/><b>Light-press:</b> Turn on/off the light intensity management function (LIM).<br/><b>Long-press:</b> Turn on/off ECO mode.</p> |
|--|--|

## 2

# Microscopy Techniques

This chapter explains the procedures for the five types of microscopy techniques supported by this microscope.

- (1) **Bright-field microscopy (See 2.1.)**
- (2) **Phase contrast microscopy (See 2.2.)**
- (3) **Simple polarizing microscopy (See 2.3.)**
- (4) **Diascopic fluorescence microscopy (See 2.4.)**
- (5) **Dark-field microscopy (See 2.5.)**

- ✓ **Turning on and off the power to the microscope**  
For detail on how to turn on and off the power to the microscope, see “3.1 Turning the Power On and Off.”

- ⓘ **Note on limits on objectives used in each microscopy technique**  
There are limits on the objectives that can be used in each microscopy technique.  
For details on the objectives and microscopy techniques, see “3.11 Combinations of Optical Parts.”

## 2.1 Bright-field Microscopy

Bright-field microscopy is the most basic observation method using diascopic illumination.

For details on the configuration and the assembly of the microscope for bright-field microscopy, see “4.2.1 Basic Configuration (For Bright-field Microscopy).”


### 1 Install the microscope.

- (1) Carry the microscope from the storage location and place it at your location.
- (2) Fix the direction of the eyepiece tube to the observation point (the binoculars are toward you).

Loosen the fixing screw for the eyepiece tube a couple of turns, and then turn the eyepiece tube.

Turn the eyepiece tube so that it faces the front, and securely tighten the fixing screw for the eyepiece tube.

- (3) Pull out the AC adapter and power cord from the storage part on the rear of the product.


Install the microscope.

### 2 Turn on the microscope.

Connect the AC adapter to power on the microscope.


- (1) Confirm that the DC plug of the AC adapter is connected to the DC input connector on the rear of the main body.
- (2) Plug the power cord of the AC adapter into an electrical outlet

Turn on the microscope.

### 3 Lower the condenser from the upper limit.


When placing the specimen, lower the condenser position by turning the condenser vertical movement knob to avoid accidental contact with the condenser lens.

Condenser vertical movement

### 4 Fully open the field diaphragm and the aperture diaphragm.

- (1) Turn the field diaphragm dial clockwise until it reaches the limit.
- (2) Set the aperture diaphragm lever of the condenser at the fully open position (leftmost).


Fully open the field diaphragm and the aperture diaphragm.

**5 Bring the 10x objective into the optical path.**

Rotate the revolving nosepiece to swing the 10x objective into the optical path.

- ✓ **Turning the nosepiece**  
The objective will click into place when rotated into position.


Set the magnification of the objective to 10x.

**6 Place the specimen on the stage, and bring the observation region into the optical path.**


- (1) Use the knob of the claw on the specimen holder to open the claw.
- (2) Place the specimen slide on the stage with the coverglass facing upward, and put the claw back to fix it.
- (3) Turn the stage knobs to bring the observation region of the specimen into the optical path.

Let illumination light go through the specimen sealed under the cover glass.

- ✓ **Using the specimen holder**  
Two glass slides can be placed in two rows on the specimen holder.

Set the specimen.

Bring the observation region into the optical path.


**7 Turn on illumination and adjust the brightness.**

The illumination on/off state and the brightness value are shown on the status display.

- (1) Press the illumination on/off button to turn on illumination.  
When illumination is turned on, the “on” symbol is shown on the status display.
- (2) Adjust the brightness by rotating the brightness control knob.

The brightness is shown in a bar graph on the status display.

- ✓ **When light intensity management is on**  
When the light intensity management function is turned on (LIM), the brightness of the selected objective is memorized. The memorized brightness is automatically loaded the next time the objective is selected.

“On” symbol      Brightness bar graph



Turn on illumination and adjust the brightness.

**8 Adjust the interpupillary distance.**

Adjust the distance between the eyepieces so that the right and left fields of view are merged into one when looking through the eyepieces.

✓ **Eye level**

When the binocular part is turned upward, the eye level is raised, and when it is turned downward, the eye level is lowered.


**Adjust the interpupillary distance.**

**9 Align the diopter ring with the standard position.**

Turn the diopter ring on the eyepieces to align its bottom edge with the engraved base line.

Perform this procedure to the right and left eyepieces.



**Set the diopter ring to the standard position.**

**10 Focus on the specimen with the 10x objective.**

Rotate the coarse and then fine focus knobs to focus on the specimen.

⚠ **Do not turn the right and left focus knobs simultaneously in opposite directions. Additionally, do not turn the coarse focus knob further after the stage has reached its upper or lower limit. This operation will damage the focusing mechanism.**

**Focus on the specimen.**


✓ **When the focal point cannot be found**

To find the focal point, rotate the stage vertical movement stopper counterclockwise to unlock the limit.


**11 Adjust the diopter.**

- (1) Swing the 40x objective in the optical path and rotate the focus knobs to focus on the specimen.
- (2) Change the objective to the 10x (or 4x) objective and focus on the specimen by rotating the diopter ring without using the focus knob.

Perform this procedure while looking into the right eyepiece with your right eye and the left eyepiece with your left eye.

(1) Focusing with the 40x objective by rotating the focus knobs

(2) Focusing with the 10x objective by rotating the diopter rings

**Adjust the diopter.**

**12 Perform focusing and centering of the condenser.**

Adjust the condenser so that illumination is correctly introduced to the specimen.


- (1) Turn the field diaphragm dial to narrow down the field diaphragm until the field diaphragm image can be seen in the field of view.
- (2) Turn the condenser vertical movement knob to focus on the field diaphragm image.

Focus on the field diaphragm image so that its outline is clearly seen.

- (3) Turn the condenser centering knob to perform centering of the condenser.


Move the condenser so that the field diaphragm image is in the center of the field of view.

- (4) Turn the field diaphragm dial to resize the field diaphragm image so that it is circumscribed with the field of view.


Narrow down the field diaphragm until the field diaphragm image can be seen.

(1) Narrow down the field diaphragm until the field diaphragm image can be seen.


(2) Focus on the field diaphragm image.  
(3) Move the field diaphragm image to the center of the field of view.

(4) Resize the field diaphragm image so that it is circumscribed with the field of view.

**Perform focusing and centering of the condenser.**


Focus on the field diaphragm image, move it to the center of the field of view, and resize the image so that it fits the field of view.

**13 Change the objective to the objective of the magnification used for observation.**

Rotate the revolving nosepiece to bring the objective of the desired magnification into the optical path.

**!** Do not directly switch from the 4x objective to the 60x or 100x objective. Be sure to focus on the specimen with the 10x objective before using the objective of a higher magnification.

- ✓ **When using an oil-immersion objective**  
Apply immersion oil between the specimen and the objective.
- ✓ **When light intensity management is on**  
When the light intensity management function is turned on (LIM), bringing the objective into the optical path loads the memorized brightness.


Switch to the objective to be used for observation.

**14 Adjust the aperture diaphragm according to the objective.**

Check the NA value of the objective, and set the aperture diaphragm lever to the position of the same figure as the NA of the objective.


→ Section 3.8, “Adjustment of Aperture Diaphragm”

Adjust the aperture diaphragm.

**15 Focus on the specimen.**


Rotate the coarse and then fine focus knobs to focus on the specimen.


Focus on the specimen.

**16 Perform focusing and centering of the condenser.**

- (1) Turn the field diaphragm dial to narrow down the field diaphragm until the field diaphragm image can be seen in the field of view.
- (2) Turn the condenser vertical movement knob to focus on the field diaphragm image.  
Focus on the field diaphragm image so that its outline is clearly seen.
- (3) Turn the condenser centering knob to perform centering of the condenser.  
Move the condenser so that the field diaphragm image is in the center of the field of view.
- (4) Turn the field diaphragm dial to resize the field diaphragm image so that it is circumscribed with the field of view.


Focus on the field diaphragm image, move it to the center of the field of view, and resize the image so that it fits the field of view.


Narrow down the field diaphragm until the field diaphragm image can be seen.

- (1) Narrow down the field diaphragm until the field diaphragm image can be seen.

- (2) Focus on the field diaphragm image.
- (3) Move the field diaphragm image to the center of the field of view.

- (4) Resize the field diaphragm image so that it is circumscribed with the field of view.

**Perform focusing and centering of the condenser.**


**17 Observe the specimen.**

Turn the stage and focus knobs to observe the specimen.

To observe the specimen with a different magnification, repeat steps 13 to 17.

**18 Turn off the microscope.**

- (1) Press the illumination on/off button to turn off illumination.  
The "on" symbol disappears from the status display.
- (2) Unplug the power cord from the electrical outlet

Turn off illumination.

**19 Store the microscope.**


After observation, turn the eyepiece tube backward and fix it.

**!** Loosen the fixing screw for the eyepiece tube a couple of turns, and then turn the eyepiece tube.

Store the AC adapter in the storage part on the rear of the microscope, wind the cords around, and then put it back in place.

**✓ AC adapter**

The AC adapter can be stored with it connected to the microscope.

**Fix the eyepiece tube backward and store the AC adapter.**

## 2.2 Phase Contrast Microscopy

Phase contrast microscopy is a method of using diasopic illumination to observe unstained, colorless and transparent specimens such as live cells.

Phase contrast microscopy uses the interference between direct light and light diffracted while passing through a colorless and transparent specimen to visualize light phase shifts as brightness variations.

Insert the phase contrast slider with the annular diaphragm into the slot on the side of the condenser, and bring the phase contrast objective with the same PH code into the optical path. The microscope is now ready for phase contrast microscopy.

For details on the configuration and the assembly of the microscope for phase contrast microscopy, see “4.2.2 Configuration for Phase Contrast Microscopy.”

### 1 Prepare the microscope.

Check that the following optical parts are usable.

#### (1) Phase contrast slider

**Slider for 10x-40x:** Supports 10x and 40x Ph objectives. When switching to phase contrast microscopy, remove the dummy slider from the slot on the side of the condenser and insert this slider.

**Slider for 100x:** Supports 100x Ph objectives. An empty position is available, so bring the empty position into the optical path in advance. When switching to phase contrast microscopy, operate the slider to bring the phase contrast position into the optical path.

#### (2) Phase contrast objective (Ph objective)

Check that the Ph objective is attached to the nosepiece.

#### ✓ Combination of optical parts




Confirm in advance that the slider is paired with a supported objective with a supported PH code.

#### (3) GIF filter (green interference filter)

Attach the GIF filter to the field lens when performing phase contrast microscopy.

#### (4) Centering telescope

Use the centering telescope when adjusting the annular diaphragm position.

### 2 Focus on the target in bright field microscopy.

Follow steps 1 to 12 of the procedure in “2.1 Bright-field Microscopy”, and focus on the target using the 10x objective with the bright-field microscopy method.

### 3 Bring the phase contrast objective into the optical path.


Turn the nosepiece to bring the phase contrast objective (Ph objective) of the desired magnification into the optical path.

The phase contrast ring inside the objective is brought into the optical path.

Check the objective magnification and the PH code.

**!** Do not directly switch from the 4x objective to the 100x objective. Be sure to focus on the specimen with the 10x objective before using the objective of a higher magnification.


- ✓ **When using an oil-immersion objective**  
Apply immersion oil between the specimen and the objective.
- ✓ **When light intensity management is on**  
When the light intensity management function is turned on (LIM), bringing the objective into the optical path loads the memorized brightness.

Switch to the phase contrast objective.

### 4 Focus on the specimen.

Rotate the coarse and then fine focus knobs to focus on the specimen.

Focus on the specimen.

### 5 Insert the phase contrast slider and bring the phase contrast position (PH) into the optical path.

#### (1) Set the aperture diaphragm lever at the fully open position.

When inserting or removing sliders, always make sure the aperture diaphragm lever is at the fully open position.


#### (2) Pull the dummy slider out of the slot on the side of the condenser.

The 100x phase contrast slider has an empty position. Bright-field microscopy can be performed by attaching the slider and bringing the empty position into the optical path in advance.

#### (3) Insert the phase contrast slider so that the centering screw of the annular diaphragm faces the front of the microscope, and the operating knob is on the left side.

Remove the slider knob, insert the slider facing the appropriate direction, and re-attach the knob to the slider.

Bring the position of the same PH code as the objective into the optical path. The annular diaphragm is brought into the optical path.

#### Inserting the phase contrast slider

- (1) **Remove the knob from the slider.**  
The knob is in the form of a screw. Unscrew the knob to remove it.
- (2) **Insert the slider into the slot on the side of the condenser.**  
The slider has protrusions to prevent incorrect insertion.
- (3) **Attach the knob to the slider you attached.**  
Attach the slider knob from the left side of the condenser.

Insert the phase contrast slider.


**6 Overwrap the annular diaphragm image with the annular image of the phase plate using the centering telescope.**

- (1) Remove one eyepiece, attach the centering telescope, and focus on the annular diaphragm image.

Focus on the specimen by turning the eyepiece while pressing down on the flange ring of the centering telescope.

- (2) Overwrap the annular diaphragm image with the annular image of the phase plate.

Insert the tool from the hole on the front of the condenser, turn the annular diaphragm centering screw of the phase contrast slider, and overlap the annular image of the phase plate of the objective with the annular diaphragm image.


Centering adjustment of the annular diaphragm


- (3) Replace the centering telescope with the eyepiece.

**7 Adjust the aperture diaphragm and field diaphragm.**


- (1) Set the aperture diaphragm lever of the condenser at the fully open position (leftmost).
- (2) Turn the field diaphragm dial to resize the field diaphragm image so that it is circumscribed with the field of view.

- (1) Focus on the annular diaphragm image using the centering telescope.

- (2) Adjust the position of the condenser so that the annular image of the phase plate and the annular diaphragm image overlap.

Adjust the field diaphragm and the aperture diaphragm.

**8 Attach the GIF filter (green interference filter) to the field lens.**

The contrast of the Ph image is improved by using the GIF filter.

**9 Observe the specimen using the phase contrast microscopy method.**

Turn the stage and focus knobs to observe the specimen.

To observe the specimen with a different magnification, repeat steps 3 to 7.

**10 Switch to bright-field microscopy.**

**(1) Remove the phase contrast slider to remove the annular diaphragm from the optical path.**

Remove the slider knob on the left side of the condenser, and pull the slider out from the right side of the condenser.

**(2) Bring the objective for bright-field microscopy into the optical path.**

## 2.3 Simple Polarizing Microscopy

Simple polarizing microscopy is an observation method for measuring polarization, refraction, and other states of the specimen by bringing the analyzer and polarizer into the microscope's optical path to limit the light's vibration direction.

For details on the configuration and the assembly of the microscope for simple polarizing microscopy, see "4.2.3 Configuration for Simple Polarizing Microscopy."

✓ **Supported polarizing microscopy**

Of the polarizing microscopy methods, this microscope only supports simple polarizing microscopy. Polarizing microscopy using sensitive color is not supported.

### 1 Prepare the microscope.

Check that the following optical parts are usable.

#### (1) Simple polarization analyzer


When switching to simple polarizing microscopy, attach it to the dovetail for the eyepiece tube mount.

(1) Simple polarization analyzer

#### (2) Simple polarization polarizer

When switching to simple polarizing microscopy, attach it to the field lens.


(2) Simple polarization polarizer

### 2 Focus on the target using the bright field microscopy method.

Follow steps 1 to 12 of the procedure in "2.1 Bright-field Microscopy", and focus on the target using the 10x objective with the bright-field microscopy method.

### 3 Bring the position without samples into the optical path.

Turn the stage knob to move the specimen and bring a portion where there is no sample under the cover glass into the optical path.

Bring the position without samples into the optical path.

#### 4 Bring the analyzer and the polarizer into the optical path.

##### (1) Turn off diascopic illumination.

Press the illumination on/off button to turn off diascopic illumination.

##### (2) Attach the simple polarization analyzer.

Remove the eyepiece tube, attach the simple polarization analyzer to the dovetail for the eyepiece tube mount, and re-attach the eyepiece tube at its original position.

Make sure the mark showing the analyzer orientation faces the same position as the screws on the rear of the microscope. The (vibration) direction of the analyzer is fixed at the microscope front and back direction.


##### (3) Attach the simple polarization polarizer.

Attach the simple polarization polarizer unit over the field lens.

Make sure the mark showing the polarizer orientation faces the front of the microscope, and keep the fixing screw in the polarizer unit loose.

##### (4) Turn on diascopic illumination.

Press the illumination on/off button to turn on diascopic illumination.

Bring the analyzer and the polarizer into the optical path.

#### 5 Adjust the angle of the polarizer so that the orientations of the analyzer and polarizer are at right angles to each other.

##### (1) Fully open the aperture diaphragm.

##### (2) Pull one eyepiece out of the eyepiece tube, and look inside the eyepiece tube sleeve.

When you look inside the eyepiece tube sleeve, you can see a black striped pattern showing the polarization state. Turning the entire polarizer unit changes the pattern.

##### (3) Turn the entire polarizer unit so that a dark cross is visible inside the eyepiece tube sleeve.

If you turn the polarizer unit so that its mark faces the front of the microscope, the polarizer is set to a horizontal orientation, and the orientations of the polarizer and analyzer converge at a right angle.


This state is known as “crossed Nicol.” If the orientations successfully converge at a right angle, a dark cross is visible inside the eyepiece tube sleeve.

##### (4) Tighten the fixing screw in the polarizer unit, and fix the polarizer in place.

The fixing screw must come in contact with a boundary protrusion of the field lens.

To fix the polarizer in an open Nicol state, remove the fixing screw, move it to another screw hole, and tighten the screw.

##### (5) Replace the centering telescope with the eyepiece.

Turn the polarizer so that a dark cross is visible.


Open Nicol state (mark on the left side)  
(To fix the polarizer in place, move the fixing screw to another position.)

**6 Bring an objective into the optical path.**

Rotate the revolving nosepiece to bring the objective of the desired magnification into the optical path.

**!** Do not directly switch from the 4x objective to the 100x objective. Be sure to focus on the specimen with the 10x objective before using the objective of a higher magnification.

- ✓ **When using an oil-immersion objective**  
Apply immersion oil between the specimen and the objective.
- ✓ **When light intensity management is on**  
When the light intensity management function is turned on (LIM), bringing the objective into the optical path loads the memorized brightness.


Switch to the desired objective.

**7 Focus on the specimen.**

Rotate the coarse and then fine focus knobs to focus on the specimen.

Adjust the brightness of diascopic illumination and move the observation region by moving the stage as required.


Focus on the specimen.

**8 Adjust the field diaphragm.**

Turn the field diaphragm dial to resize the field diaphragm image so that it is circumscribed with the field of view.

- ✓ **Adjustment of field diaphragm**  
Each time you change the objective, adjust the field diaphragm.

Adjust it so that the image is circumscribed with the field of view

Adjust the field diaphragm.

**9 Observe the specimen using the simple polarizing microscopy method.**

Turn the stage and focus knobs to observe the specimen.

To observe the specimen with a different magnification, repeat steps 6 to 8.

**10 Switch to bright-field microscopy.**

**(1) Turn off diascope illumination.**

Press the illumination on/off button to turn off diascope illumination.

**(2) Remove the simple polarization analyzer.**

Remove the eyepiece tube, remove the simple polarization analyzer from the dovetail for the eyepiece tube mount, and re-attach the eyepiece tube at its original position.

**(3) Remove the simple polarization polarizer.**

Loosen the fixing screw in the polarizer unit, and remove the polarizer unit from the field lens.

**(4) To continue observation, turn on diascope illumination.**

Press the illumination on/off button to turn on diascope illumination.

## 2.4 Diascopic Fluorescence Microscopy

Diascopic fluorescence microscopy is a microscopy method for using diascopic illumination to examine specimens stained with fluorescent dye or fluorescent protein. You can observe excited fluorescence as an image by directing diascopic illumination of a particular wavelength at objects stained with special fluorescent dye.

This microscope can be used for observation by diascopic fluorescence (GFP-B) by using the E2-F-FL Dia-FL Set GFP-B. You can perform observation using GFP, FITC, Alexa 488, and more.

✓ **Note on the brightness of excited fluorescence**

The brightness of fluorescent light excited during diascopic fluorescence microscopy might be reduced, for example due to the expression efficiency of the sample GFP.

Nikon recommends darkening the room because images are not clearly visible if the ambient light is too bright.

For details on the configuration and the assembly of the microscope for diascopic fluorescence microscopy, see “4.2.4 Configuration for Diascopic Fluorescence Microscopy.”

### 1 Prepare the microscope.

Attach the following optical parts and the shielding plate to the microscope in advance:

#### (1) EX filter slider

When performing bright-field microscopy, bring the BF position into the optical path, and when performing diascopic fluorescence microscopy, bring the GFP position into the optical path.

#### (2) BA filter slider


When performing bright-field microscopy, pull the slider out to the click position and bring the empty position into the optical path. When performing diascopic fluorescence microscopy, push the slider in all the way and bring the BA filter into the optical path.

#### (3) Diascopic fluorescence objective

When performing diascopic fluorescence microscopy, you can use three types of objectives: CFI E Plan 4x, CFI E Plan 10x, and CFI E Plan 40x.

#### (4) Shielding plate

Attach the shielding plate to prevent strong light including UV light from the specimen through the condenser from flashing into your eyes.

### 2 Focus on the target using the bright field microscopy method.


Follow steps 1 to 12 of the procedure in “2.1 Bright-field Microscopy”, and focus on the target using the 10x objective with the bright-field microscopy method.

**3 Bring the diascope fluorescence objective into the optical path.**

Turn the nosepiece to bring the objective into the optical path.

- ✓ **When light intensity management is on**  
When the light intensity management function is turned on (LIM), bringing the objective into the optical path loads the memorized brightness.

- ⚠ **Note on using light intensity management**  
If diascope fluorescence microscopy was being performed last time the microscope was used, beware that bright light might enter the optical path.


Switch to the diascope fluorescence objective.

**4 Focus on the specimen.**

Rotate the coarse and then fine focus knobs to focus on the specimen.

Adjust the brightness of diascope illumination and move the observation region by moving the stage as required.


Focus on the specimen.

**5 Adjust the field diaphragm.**

Turn the field diaphragm dial to resize the field diaphragm image so that it is circumscribed with the field of view.

- ✓ **Adjustment of field diaphragm**  
Each time you change the objective, adjust the field diaphragm.

Adjust it so that the image is circumscribed with the field of view

Adjust the field diaphragm.

**6 Bring the excitation filter and barrier filter into the optical path, and adjust the illumination to set the microscope to diascope fluorescence microscopy.**

- (1) **From the right side, push the EX filter slider inserted in the slot on the side of the condenser. This sets the slider to the GFP position.**

The excitation filter is brought into the optical path.

- !** **Move your eyes away the eyepieces during position switching.**  
**Fluorescence microscopy is performed with diascope illumination set to maximum. Therefore, bright light might reach the binoculars when the position of the fluorescence microscopy slider is switched.**  
**Move your eyes away from the binoculars during position switching.**

Push the slider into the click position. The aperture diaphragm lever is automatically set to fully open (NA 1.25), and is no longer operable.

- (2) **Push the BA filter slider all the way into the front slot of the nosepiece.**

The barrier filter is brought into the optical path.

- (3) **Adjust the brightness of diascope illumination.**

Brighten the diascope illumination so that the specimen is excited and emits fluorescence.


- ✓ **Illumination during diascope fluorescence microscopy**

To perform fluorescence microscopy, use illumination close to the maximum illumination intensity.


**7 Observe the specimen using the diascope fluorescence microscopy method.**

Turn the stage and focus knobs to observe the specimen.

To observe the specimen with a different magnification, repeat steps 6 to 8.

- (1) Set the EX filter slider to the GFP position.  
 (2) Push the BA filter slider in all the way.

- (3) Adjust the brightness of diascope illumination.

**Set the microscope to use diascope fluorescence.**

**8 Switch to bright-field microscopy.**

**(1) Darken diascope illumination, and then turn off diascope illumination.**

Turn the brightness control knob to darken diascope illumination. Then, press the illumination on/off button to turn off diascope illumination.

**(2) Pull the BA filter slider up to the click position.**

The barrier filter is removed from the optical path.

**!** Move your eyes away the eyepieces during position switching.  
Beware that if you remove the barrier filter from the optical path while looking into the eyepiece during diascope fluorescence microscopy, light might flash into your eye, causing eye damage.  
When setting the microscope back to bright-field microscopy, either turn off diascope illumination or keep your eyes away from the eyepiece.

**(3) Insert the EX filter slider from the left side up to the click position, and set it at the BF position (for bright-field microscopy).**

The empty position is brought into the optical path.

**(4) To continue observation, turn on diascope illumination.**

Press the illumination on/off button to turn on diascope illumination.

## 2.5 Dark-field Microscopy

Dark-field microscopy is a method of using oblique illumination to observe unstained specimens. A bright observation target can be visualized on a dark background in high contrast.

In dark-field microscopy, the annular diaphragm narrows the illumination to an annular shape, and light passes through the edge of the condenser lens to obliquely illuminate the specimen. Using an objective with a smaller numerical aperture than the condenser lens allows only light scattered by the specimen to pass through the objective to form a dark-field image.

For details on the configuration and the assembly of the microscope for dark-field microscopy, see “4.2.5 Configuration for Dark-field Microscopy.”

### 1 Prepare the microscope.

Check that the following optical parts are usable.

#### (1) Slider for dark-field microscopy


When performing bright-field microscopy, bring the empty position into the optical path, and when performing dark-field microscopy, bring the dark-field position into the optical path.

#### (2) Objective for dark-field microscopy

Check that the objective is attached to the nosepiece.

#### ✓ Combination of optical parts

There are limits on the NA values of objectives that can be used with dark-field microscopy.


### 2 Focus on the target using the bright field microscopy method.

Follow steps 1 to 12 of the procedure in “2.1 Bright-field Microscopy”, and focus on the target using the 10x objective with the bright-field microscopy method.

### 3 Insert the dark-field slider and bring dark-field position into the optical path.

Insert the dark-field slider from the right side up to the click position, and bring the dark-field position into the optical path.

The dark-field annular diaphragm is brought into the optical path. The microscope is now ready for dark-field microscopy.


Bring the objective for dark-field microscopy into the optical path.

**4 Fully open the field diaphragm and the aperture diaphragm.**

- (1) Turn the field diaphragm dial clockwise until it reaches the limit.
- (2) Set the aperture diaphragm lever of the condenser at the fully open position (leftmost).

✓ **Condenser position**

Move the condenser to the position where the image is clear under dark-field illumination. The position slightly lower than the upper limit would be suitable.

**Fully open the field diaphragm and the aperture diaphragm.**

**5 Bring an objective into the optical path.**

Rotate the revolving nosepiece to bring the objective of the desired magnification into the optical path.


❗ **Do not directly switch from the 4x objective to the 100x objective. Be sure to focus on the specimen with the 10x objective before using the objective of a higher magnification.**

✓ **When using an oil-immersion objective**

Apply immersion oil between the specimen and the objective.

✓ **When light intensity management is on**

When the light intensity management function is turned on (LIM), bringing the objective into the optical path loads the memorized brightness.


**Switch to the suitable objective.**

**6 Adjust the brightness of diascope illumination.**

Adjust the brightness of diascope illumination by turning the brightness control knob.

In dark-field microscopy, the background darkens, and the target is illuminated from an angle. Adjust the brightness to observe the specimen as desired.

Brightness control knob

**Adjust the brightness of diascope illumination.**


**7 Focus on the specimen.**

Rotate the coarse and then fine focus knobs to focus on the specimen.

**8 Observe the specimen using the dark-field microscopy method.**

Turn the stage and focus knobs to observe the specimen.

To observe the specimen with a different magnification, repeat steps 5 to 7.

**Focus on the specimen.**

**9 Switch to bright-field microscopy.**

**(1) Darken diascope illumination, and then turn off diascope illumination.**

Turn the brightness control knob to darken diascope illumination. Then, press the illumination on/off button to turn off diascope illumination.

**(2) Insert the dark-field slider from the left side up to the click position, and bring the empty position into the optical path.**

The annular diaphragm is removed from the optical path, and the microscope is set to bright-field microscopy.

**!** Move your eyes away the eyepieces during position switching.  
Beware that if you remove the barrier filter from the optical path while looking into the eyepiece during dark-field microscopy, bright light might enter your eyes.  
When setting the microscope back to bright-field microscopy, either turn off diascope illumination or keep your eyes away from the eyepiece.

**(3) To continue observation, turn on diascope illumination.**

Press the illumination on/off button to turn on diascope illumination.

## 3

## Detailed Explanation

## 3.1 Turning the Power On and Off

## 3.1.1 Turn On/Off the Power to the Microscope

**Powering on the microscope**

Connect the AC adapter to turn on this microscope.

Connect the DC plug of the AC adapter to the DC IN connector on the rear of the product. Plug a power cord into the AC adapter's inlet for equipment, and the wall outlet.

When the power is turned on, the startup screen appears on the status display. Then, the control screen appears as usual.

**Powering off the microscope.**

To turn off the microscope, unplug the AC adapter.

Startup screen

Control screen

## 3.1.2 ECO Mode and SLEEP Mode

The ECO and SLEEP modes are used to put the microscope into the energy-saving state. During the energy-saving state in both modes, diascopic illumination and the LCD screen are turned off to save energy.

**ECO mode:**


When the microscope (turned on) is left idle for a specified period of time, the microscope is put into the energy-saving state. Long-press the brightness control knob to turn on or off the ECO mode (default: on). When the ECO mode is turned on, "ECO" is shown on the status display.

✓ **Setting the time to put the microscope into the energy-saving state.**

You can set how long the ECO mode waits before the microscope is put into the energy saving state from the SET UP screen (default: 15 minutes).

✓ **Observation through only eyepieces or long-term observation (without part operation)**

Turn off the ECO mode.

**Operating the brightness control knob**

- **Long-press:** Turn the ECO mode on or off.

**Operating the illumination on/off button**

- **Long-press:** Activate the SLEEP mode.

**SLEEP mode:**


The microscope can be put into the energy saving state anytime.

Long-press the illumination on/off button to immediately put the microscope into the energy-saving state.

**Recovery from the energy-saving state**

Perform the following to bring the microscope back to the normal operating state:

- Turn or press the brightness control knob.
- Press the illumination on/off button.
- Turn the nosepiece.

Control screen

## 3.2 Operating Diascopic Illumination

This microscope contains a diascopic illumination device using an LED light source.

### 3.2.1 Turning Diascopic Illumination On and Off and Adjusting the Brightness

#### Turning diascopic illumination on and off


To turn diascopic illumination on or off, use the illumination on/off button on the left side of the microscope.

Pressing the button once toggles illumination on or off. When illumination is on, the “on” symbol is shown on the status display.

#### Adjusting the brightness of diascopic illumination (brightness adjustment)

To adjust the brightness of diascopic illumination (brightness adjustment), turn the brightness control knob.

The illumination brightness (brightness state) is shown in a bar graph on the status display.


#### Operating the illumination on/off button

- **When pressed:** Toggles diascopic illumination on or off

#### Operating the brightness control knob

- **When turned:** Adjusts the brightness of diascopic illumination


### 3.2.2 Using Light Intensity Management

This microscope features a light intensity management function (LIM).

To turn light intensity management on or off, briefly press the brightness control knob (the default setting is “on”).


When the light intensity management function is on, the brightness value for each position in the nosepiece is saved to the microscope's memory. The next time a position is brought into the optical path, the brightness value from the last time that position was used is automatically loaded.

When light intensity management is on, the “LIM” symbol is shown on the status display on the front of the microscope.

#### Operating the brightness control knob

- **When briefly pressed:** Turns light intensity management (LIM) on or off


### 3.3 Operating the Stage Specimen Holder

#### 3.3.1 Placing a Specimen Using the Specimen Holder

Specimens are placed using the specimen holder installed on the stage.

- (1) Use the knob of the claw on the specimen holder to open the claw.
- (2) Place the specimen slide on the stage with the cover glass facing up.
- (3) Fix the specimen in place by returning the claw to its original position.

- ✓ **Placing a specimen**  
Two glass slides can be mounted on the specimen holder.
- ✓ **Cleaning the specimen**  
Clean the surface of the glass slide and the cover glass as necessary.

Placing a specimen using the specimen holder


#### 3.3.2 Operating the Stage

You can set the specimen fixed to the specimen holder to the desired position by operating the stage motion knobs to move the stage back and forth, and right and left. Set the specimen so that the region you want to observe is in the optical path.


The stage is capable of the following travel strokes:

- **Forward and back (top knob):** 52 mm
- **Left and right (bottom knob):** 76 mm
- ✓ **Travel scale**  
The stage has the scales that indicate the travel of the specimen in the horizontal and vertical directions. You can use the main scale and vernier scale to read the travel in units of 0.1 mm.

The end of the travel scale shows the icon of the stage motion knob to be operated.

Operating the stage


## 3.4 Operating the Nosepiece and Objectives

### 3.4.1 Switching Objectives

To switch objectives to use for observation, manually turn the turret ring of the nosepiece right or left until it reaches the click position.

Objectives are usually set in such an order that the objective magnification increases as the nosepiece is turned clockwise (as viewed from the top of the nosepiece).

- ✓ **When light intensity management is on**  
When the light intensity management function is turned on (LIM), bringing the objective into the optical path loads the memorized brightness.

Operating the nosepiece

### 3.4.2 Oil-Immersion Observation

The “Oil” mark on the side of an objective indicates that it is an oil-immersion type objective. (The oil-immersion objective also has a black band around the barrel end.) An oil-immersion objective is used with the immersion oil applied between the front of the objective and the coverglass.

For an oil-immersion objective with a numerical aperture of 1.0 or more, also apply immersion oil to the condenser to take full advantage of its performance. (An oil-immersion type condenser needs immersion oil to be applied between the front of the condenser and the specimen.)

The condenser has an oil receptacle around the front of its lens.

#### Example of oil-immersion

Apply immersion oil to the condenser and then objective.

- **Condenser:**  
Add a drop of oil on the condenser lens, and then place the specimen on the stage.
- **Objective:**  
Rotate the revolving nosepiece to move the objective out of position. Add a drop of oil to the specimen. Slowly rotate the revolving nosepiece to bring the objective back into position.

#### Eliminating air bubbles

Make sure that air bubbles are not trapped during oil application. Air bubbles degrade the image.

Do any of the following to eliminate air bubbles:

- Rotate the revolving nosepiece to move the objective back and forth once or twice.
- Add another drop of oil.
- Wipe off the oil and apply again.

## Handling of the immersion oil

- Use a minimum quantity of oil. If too much oil is applied, surplus oil may flow out onto the stage and the condenser and degrade performance.
- Any oil residue left on the lenses of oil-immersion type objectives or adhesion of oils on the front lens of dry type objectives will degrade image quality. After completing oil-immersion observation, be sure to clean the objective, condenser, and any other parts that may be stained by oil.
- To wipe surplus oil on the front of the condenser lens, remove the condenser.
- Use petroleum benzine to wipe off oil and finish with absolute alcohol (ethyl or methyl alcohol). If petroleum benzine is not available, use methyl alcohol instead. In that case, wipe off the oil several times as the detergency of methyl alcohol is weaker than petroleum benzine. (You can wipe off the oil generally 3 to 4 times.)
- When handling petroleum benzine and absolute alcohol, be sure to follow the instructions described in “4. Handling flammable solvents” in “WARNING” at the beginning of this manual.
- Close the oil container cap tightly after use. Make sure that the cap is closed tight after refilling the container. Check the cap periodically to make sure it has not come loose, allowing oil to leak out.
- Do not press the container hard. Oil may splash out.
- If you find an oil drips around the container, wipe them off.
- Avoid contact of immersion oil with eyes or skin. In the event of contact with eyes or skin, take one of the following measures although Nikon immersion oil does not contain any toxic ingredients.
 

<b>Contact with skin:</b>	<b>Rinse your skin thoroughly with soap and water.</b>
<b>Contact with eye:</b>	<b>Rinse your eye thoroughly with water (more than 15 minutes) and see a doctor.</b>
- Do not leave immersion oil in the sun (ultraviolet rays may damage it).

### 3.4.3 Performing Water-Immersion Microscopy

The “WI” or “W” mark on the side of an objective indicates that it is a water-immersion-type objective. When using such an objective, perform water immersion (using distilled water or physiological saline) between the specimen and the front of the objective. Operate the microscope in the same way as for oil immersion.

Because water evaporates easily, make sure that the water is still present during observation. Beware that using large amounts of water might cause the water to overflow into the stage or condenser, causing corrosion.

#### ✓ **Cleaning up moisture**

After using the microscope, wipe any moisture off the surface of the objective. Then, perform a finishing wipe with absolute alcohol.


If evaporated moisture leaves a stain, gently wipe off the stain using a small amount of neutral detergent. Then, finish cleaning the surface by performing a finishing wipe with absolute alcohol.

## 3.5 Focusing Operations

### 3.5.1 Focusing on a Specimen

This procedure is to operate the fine and coarse focus knobs to raise and lower the stage so that the specimen to be observed can be seen most clearly.

The UP and DOWN indications on the sides of the microscope base indicate the direction in which the stage moves by rotation of the focus knobs. For both the fine and coarse focus knobs, turning in the UP direction raises the stage, and turning in the DOWN direction lowers the stage.

#### ✓ Fine scale

One rotation of the fine focus knob moves the stage vertically by 0.2 mm. Fine tick marks are placed at 2- $\mu$ m intervals of stage travel.

### Standard procedure

- (1) Swing the 10x (or 4x) objective into the optical path.
- (2) Rotate the coarse focus knob to raise the stage to its upper limit.

Since the working distances of 4x and 10x objectives are long, these objectives do not touch the specimen even when the stage is raised to its upper limit provided that the slide and coverglasses of a standard thickness are used. (The standard thickness for slide is 1.2 mm and that for coverglass is 0.17 mm.)


- (3) While looking into the eyepiece, slowly turn the coarse focus knob to lower the stage. When you find a position where the specimen is roughly in focus, let go of the coarse focus knob.

⚠ When rotating the coarse focus knob while looking into the eyepiece, be sure to rotate it in the direction that lowers the stage.

#### ✓ When the focal point cannot be found

To find the focal point, rotate the stage vertical movement stopper counterclockwise to unlock the limit.

- (4) Rotate the fine focus knob and precisely focus on the image.

## Using the objective of a higher magnification for observation

First, use the 10x objective to focus on the specimen, and then change it to the desired objective and rotate the fine focus knob to precisely focus on the specimen.


**!** When operating the focus knob while using the objective of a higher magnification, care must be taken to prevent the slide glass from touching the objective. When you raise the stage using the coarse focus knob, look at the gap between the upper surface of the specimen and the front of the objective from the side.

## Using working distances

The stage can be raised by using the working distance of the objective.

While looking at the microscope from the side, carefully rotate the coarse focus knob to bring the specimen close to the objective. When the distance between the specimen and the front of the objective becomes slightly shorter than the working distance, let go of the coarse focus knob. The specimen is now almost in focus.

Look into the eyepieces and rotate the fine focus knob in the direction that lowers the stage to precisely focus on the specimen.

## Focusing by switching the objective from 40x to 10x

Swing the 40x objective into the optical path. While looking at the microscope from the side, rotate the coarse focus knob until the specimen almost touches the objective (about 0.5 mm apart from the front of the objective).


Switch to the 10x objective, look into the eyepieces, and rotate the fine focus knob slightly to find the focal point. Be careful not to hit the objective with the specimen.

### 3.5.2 Torque Adjustment of the Coarse Focus Knob

The tension (torque) of the coarse focus knob rotation can be adjusted.

To increase the tension, turn the coarse focus knob torque adjustment ring counterclockwise.

To decrease the tension, turn the ring clockwise. Do not decrease the tension too much. If it is too loose, the stage will fall under its own weight.

### 3.5.3 Replacing a Specimen Using the Stage Vertical Movement Stopper

When focusing on the specimen with 40x or higher magnification objective, you will find the specimen is brought very near to the objective (\*1).

At this time, special care must be taken to avoid damaging the specimen when replacing it. In a case like this, use the stage vertical movement stopper for easy specimen replacement.

The stage vertical movement stopper has the labels "FOCUS LIMIT", "HIGH", and "LOW", which show the direction of rotation and vertical movement. Operate the stage vertical movement stopper according to the following steps.




- (1) Turn the stage vertical movement stopper counterclockwise (in the HIGH direction) until it stops.
- (2) Focus on the specimen, and then turn the stage vertical movement stopper clockwise (in the LOW direction) until it stops. Rotating the stage vertical movement stopper too much causes the set focal point to shift.
- (3) Performing the above steps sets the stage so that it does not move above that position.
- (4) Lower the stage by the focus knob, and replace the specimen. After that, when turning the focus knob, you can find that the stage stops near the position set in (2).
- (5) If the specimen is out of focus, turn the stage vertical movement stopper counterclockwise (in the HIGH direction) to release the set position, and then slightly operate the fine focus knob to adjust the focus.

\*1: The distance between the front of the objective and the specimen when the specimen is in focus is called the "working distance" of the objective. See Section 3.10, "Microscope Terminology."


#### ✓ Operating the stage vertical movement stopper

Turn the stopper behind the microscope arm by fingers as shown in the figure below on the left.



If you cannot reach the stopper in this way, turn the stopper from the rear of the microscope arm as shown in the figure below on the right.

Turn the stopper until it stops (in the HIGH direction)

Focusing

Turn the stopper until it stops (in the LOW direction)

When raising the stage by the focus knob, the stage stops near the focused position. The stage does not move above that position.

## 3.6 Operating the Binoculars

### 3.6.1 Adjustment of Diopter


This procedure is to adjust the diopter ring on the eyepieces according to the difference between your right and left eyesight.

✓ **Effect of diopter correction**


This adjustment enables the user to take full advantage of objective performance, and reduce loss of focus when switching magnifications.

- (1) Turn each diopter ring on the right and left eyepieces to align the diopter ring's bottom edge with the engraved base line.

This position is the base position for diopter adjustment.


- (2) Swing the 40x objective in the optical path. Rotate the coarse and then fine focus knobs to bring the specimen in focus. See the description in "■ Focusing by switching the objective from 40x to 10x" in Section 3.3, "Focusing."


- (3) Swing the 10x (or 4x) objective into the optical path.

- (4) While looking into the right eyepiece with your right eye, focus on the specimen by rotating the right diopter ring without using the focus knob.

- (5) While looking into the left eyepiece with your left eye, focus on the specimen by rotating the left diopter ring without using the focus knob.
- (6) Repeat steps (2) to (5) to check the specimen is in focus.

### 3.6.2 Adjusting Interpupillary Distance

“Adjusting interpupillary distance” means to adjust the distance between the eyepieces to match the distance between the observer’s eyes. Performing this adjustment makes it easier to observe specimens with both eyes.

#### Selecting an eye point



The binoculars of this microscope can be turned both upward and downward. To raise the eye point, turn the binoculars upward. To lower the eye point, turn the binoculars downward.

#### Adjusting interpupillary distance

After completing diopter adjustment of the eyepieces, bring the 10x objective into the optical path and focus on the specimen. Then, adjust the interpupillary distance of the binoculars so that looking into the eyepieces with both eyes results in a unified field of view.

#### ✓ Using the interpupillary distance scale lines

The binoculars are labeled with interpupillary scale lines. If you remember your own interpupillary distance, you can easily adjust the interpupillary distance the next time you use the microscope.

Adjusting interpupillary distance

## 3.7 Operating the Condenser


### 3.7.1 Focusing the Condenser

Adjust the vertical position of the condenser so that the field diaphragm image is clearly visible in the field of view.

**!** **Caution: Do not raise the condenser too much.**

If you raise the condenser too much, the condenser lens might hit the specimen's glass slide, causing damage to the condenser lens or the glass slide.

When raising or lowering the condenser lens, be careful not to hit the specimen's glass slide.

### 3.7.2 Centering the Condenser


Adjust the position of the condenser so that the illumination is correctly centered in the field of view. This operation is called “centering the condenser.”

You can adjust the position of the condenser by turning the condenser centering knobs on each side of the aperture diaphragm lever.

#### Adjusting the position of the condenser

Bring the 10x objective into the optical path, narrow the field of view, and adjust the position of the condenser so that the field diaphragm image is centered in the field of view.

When performing strict position adjustment, bring a higher-magnification objective into the optical path and perform centering as usual.


### 3.8 Adjustment of Aperture Diaphragm

The size of the condenser aperture diaphragm must be adjusted according to the NA value of the objective in the optical path.

The condenser is labeled with an adjustment position mark for the aperture diaphragm, which serves as a guideline for the adjustment position. Turn the aperture diaphragm lever to the same number as the NA value of the objective in the optical path.

Each time you change the objective magnifications, align the aperture diaphragm lever to the same figure as the magnifying power of the objective.

Adjustment position mark for the aperture diaphragm

#### ✓ Aperture Diaphragm

- The aperture size is increased or decreased by turning the condenser aperture diaphragm lever. If the aperture diaphragm is closed, the brightness and resolution are decreased but the contrast and range of focus are increased. If the aperture diaphragm is opened, the brightness and resolution are increased but the contrast and range of focus are decreased. Since the image resolution will be degraded when the aperture diaphragm is closed too much, do not close the aperture diaphragm too much except when observing a specimen with extremely low contrast, such as a near-transparent specimen.
- The aperture diaphragm controls the numerical aperture of the illumination. Do not use it to control brightness. Use the brightness control knob to control the brightness.

- E PLAN → Objective type
- 40x/0.65 → Magnification 40 times/Numerical aperture 0.65
- $\infty$ /0.17 → Mechanical tube length  $\infty$ /Coverglass thickness 0.17 mm
- WD 0.65 → Working distance (Distance between the tip of the objective and the specimen) 0.65 mm

#### Specifications of objective


### 3.9 Adjustment of field diaphragm

“Adjusting the field diaphragm” means narrowing the illumination to limit the field of observation of the specimen.

To resize the field diaphragm, turn the field diaphragm dial. Usually, observation is performed with the field diaphragm narrowed until it roughly circumscribes the edge of the field of view.

Note that the field diaphragm cannot be narrowed to inside field of view when using a 60x or higher magnification objective.

- ✓ **Optimum size of the field diaphragm**  
Usually, the optimum size is when the field diaphragm roughly circumscribes the edge of the field of view.  
If the field diaphragm is opened too wide, causing an unnecessarily wide field of illumination, stray light enters from the outside, resulting in flares and reducing the contrast of the optical image.  
This also causes the specimen's discoloration field to become too wide.
- ✓ **Adjustment of field diaphragm**  
Each time you change the objective, adjust the field diaphragm.

**Adjustment of field diaphragm**

## 3.10 Tips for Each Microscopy Method

### 3.10.1 Tips for Phase Contrast Microscopy

Phase contrast microscopy is suited for observing colorless and transparent specimens, unstained or lightly stained specimens, specimens that have lost their color, and ultra-thin sections for electron microscopes. Shaded or strongly stained specimens are not suited to phase contrast microscopy.

The appearance of the phase contrast image is influenced by factors such as the phase difference or shape of the specimen and the characteristics of the objective. Be careful of the following when preparing specimens and selecting Ph objectives.

#### Select specimens that do not de-center the Ph annular diaphragm.

When light-scattering specimens or specimens that have a lens or prismatic effect are observed, the Ph annular ring becomes de-centered. Special care is needed with specimens such as thick live specimens, coarse specimens, and specimens using microplates, which affect the lens or prism and de-center the Ph annular ring, worsening the image.

#### Ph objectives and specimens

Ph objectives can be divided into “achromat”, “plan achromat”, “plan fluor”, and “plan apochromat” objectives depending on their chromatic aberration and degree of field curvature correction. These lenses can also be divided into several types depending on the characteristics of their internal phase ring. Satisfactory results cannot be achieved if the specimen's degree of phase difference or the characteristics of the phase ring are unsuitable. For details on the characteristics of the use of phase contrast objectives, see the following table.

- When using a Ph objective for dark contrast, make sure the specimen's phase difference does not exceed the objective's degree of tolerance (latitude) for phase difference. If the specimen's degree of phase difference exceeds the degree of tolerance for phase difference, the image shines more brightly than the background, making observation impossible.
- When preparing a phase-difference specimen, you can adjust the phase difference by changing the refractive index, for example by changing the specimen thickness, mounting agent, or culture fluid.
- Even specimens that appear in low contrast when observed using a DLL objective might yield good results when observed using a DM objective.

Characteristics of the use of phase contrast objectives

Phase contrast objective		Appearance of image	Contrast		Latitude	Example of application
Dark contrast	DLL DL	Generally, objects with a large phase difference appear dark. Accordingly, the image appears blackened in a comparatively bright field of view, resulting in a similar image to bright-field microscopy.	Dark contrast is suited to detailed observation, mainly via micro-contrast.	Halftones (wide scope of use)	Low- and mid-range phase difference and absorptive objects (stained objects)	Fungal spores, general live cells, slightly thick specimens, bacteria, stained specimens, insect eggs, fat globules, crystals, and the like
	DM			High-contrast tones (comparatively narrow scope of use)	Low-range transparent objects	Bacterial or protozoan flagella, fibrin fibers, fine granules, sections with appropriately selected mounting agents, ultra-thin sections, and the like
Bright contrast	BM	Generally, objects with a large phase difference appear bright. Accordingly, the image appears brightened in a comparatively dark field of view, resulting in a similar image to dark-field microscopy.	Bright contrast is suited to observing the forms of detailed fibers or granules, detection, or calculation, mainly via macro-contrast.		Mostly full-range	Bacterial or protozoan flagella, fibrin fibers, fine granules, blood-count calculation, and the like

## Notes on limits on objectives

---

When a phase contrast slider's bright-field position is used, the following limits apply depending on the objective:

- For objectives of NA 0.75 or higher, part of the field of view is lost (vignetting occurs).
- For objectives of NA 0.95 or higher, loss of field of view caused by vignetting, or reductions in peripheral light intensity occur.

Note that when a phase contrast slider's bright-field microscopy position is used, the following limits apply depending on the objective:

- For PH sliders (10x/40x): PH1: NA 0.25 - 0.5, PH2: NA 0.55 - 0.85
- For PH sliders (100x): PH3: NA 0.9 - 1.4

## Phase contrast slider

---

Phase contrast microscopy requires a phase contrast slider with a mounted Ph annular diaphragm. A phase contrast result can be obtained by correctly matching the Ph annular diaphragm and the phase ring of the Ph objective.

### Ph code

Ph objectives are labeled with a Ph code according to the size of the phase ring. The code can be [Ph1], [Ph2], or [Ph3]. To use a Ph objective, bring the Ph annular diaphragm that has the same code as the Ph objective's Ph code into the optical path. If you use a Ph annular diaphragm with a different Ph code from the Ph objective, a phase contrast result cannot be obtained.

### Aligning an annular diaphragm and phase ring

The phase contrast image is subtly affected by overlap between the annular diaphragm and phase ring. When performing strict observation or photography, confirm that the annular diaphragm and phase ring are concentric for each magnification.

Adding slight eccentricity to the annular diaphragm and phase ring creates a shadowing effect, producing a three-dimensional image. Use this technique as required by the specimen.

## Using GIF filters

---

The contrast of a phase contrast image is improved by bringing a GIF filter (green interference filter) into the optical path. GIF filters are installed on top of the field lens.

### 3.10.2 Tips for Simple Polarizing Microscopy

Simple polarizing microscopy is an observation method for measuring polarization, refraction, and other states of the specimen by bringing the analyzer and polarizer into the microscope's optical path to limit the light's vibration direction.

Of the polarizing microscopy methods, this microscope only supports simple polarizing microscopy. Polarizing microscopy using sensitive color is not supported.

## Operating the polarizer

---

- Setting the polarizer to a horizontal orientation causes the analyzer and polarizer orientations to converge at a right angle. This state is known as "crossed Nicol."
- When the polarizer is removed and only the analyzer is brought into the optical path, this state is known as "open Nicol."

### 3.10.3 Tips for Diascopic Fluorescence Microscopy

Diascopic fluorescence microscopy is a microscopy method for using diascopic illumination to examine specimens stained with fluorescent dye or fluorescent protein. You can observe excited fluorescence as an image by directing diascopic illumination of a particular wavelength at objects stained with special fluorescent dye.

#### Purpose of diascopic fluorescence microscopy

This microscope can be used for observation by diascopic fluorescence (GFP-B) by using the E2-F-FI diascopic fluorescence filter set (GFP-B). You can perform observation using GFP, FITC, Alexa 488, and more.

✓ **Note on the brightness of excited fluorescence**

The brightness of fluorescent light excited during diascopic fluorescence microscopy might be reduced, for example due to the expression efficiency of the sample GFP.

Generally, excited fluorescent light is extremely faint.

When performing fluorescence microscopy, darken the area around the microscope, for example by turning down the ambient light so that you can observe the specimen correctly.

#### Note on eye protection during slider operation

Move your eyes away from the binoculars during position switching of EX filter sliders.

Fluorescence microscopy is performed with diascopic illumination set to maximum. Therefore, bright light might reach the binoculars when the position of the fluorescence microscopy slider is switched.

**Do not remove the slider.**

**Never remove the EX filter slider from the microscope during diascopic fluorescence microscopy.**

#### Preventing UV light (using a shielding plate)

Use the shielding plate to prevent strong light including UV light from the specimen through the condenser from entering your eyes.

#### Protecting specimens and preventing discoloration

Brightly illuminating a specimen for fluorescence microscopy for an extended time might cause damage or discoloration of the specimen.

When stopping diascopic fluorescence microscopy, either turn off the illumination or lower the brightness. Make this a habitual practice.

#### Using non-fluorescent glass slide, a cover glass, and immersion oil

To obtain an image with good contrast during fluorescence microscopy, make sure that you use non-fluorescent glass slide, a cover glass, and an immersion oil specified by our company.

#### Limiting illumination to the field of observation (adjusting the field diaphragm)

“Adjusting the field diaphragm” means narrowing the illumination to limit the field of observation of the specimen.

Usually, observation is performed with the field diaphragm narrowed until it roughly circumscribes or inscribes the edge of the field of view. If the field diaphragm is opened too wide, causing an unnecessarily wide field of illumination, stray light enters from the outside, resulting in flares and reducing the contrast of the optical image. This also causes the specimen's discoloration field to become too wide.

Each time you change the objective, adjust the field diaphragm.

#### Finding targets inside the specimen

For fluorescence microscopy, the general practice is to first use phase difference to find a target and then switch to using fluorescence microscopy. However, when finding targets using bright-field microscopy, you will need to take the following kinds of measures:

- In bright-field microscopy, first use a 10x objective, and then narrow the condenser as required.
- If raising the magnification makes it difficult to find the target, switch to fluorescence microscopy, and then find the target while suppressing the excitation light.

Another method is to work while looking at the edge plane of the cover glass to identify the rough position of the target.

### 3.10.4 Tips for Dark-field Microscopy

#### Notes on limits on objectives

---

Objectives that can be used for dark-field microscopy are limited to the following:

- **For dark-field microscopy:** NA 0.65 or lower
- **For observation using a bright-field position:** NA 1.25 or lower

✓ **Note on the bright-field position of the dark-field slider**

Objectives that can be used in the bright-field position of the dark-field slider have different conditions than for regular bright-field microscopy.

With objectives of NA 1.0 or higher, vignetting might occur in the periphery, causing the image to darken.

#### Note on eye protection during slider operation

---

Move your eyes away from the binoculars during position switching of dark-field sliders.

When you switch from a dark field to a bright field, bright light might reach the binoculars.

## 3.11 Combinations of Optical Parts

### Compatibility between objectives and microscopy methods

#### Compatibility between optical parts and microscopy methods

Objective	Magnification	Bright-field	Phase Contrast	Simple Polarizing	Diascopic Fluorescence	Dark-field
E Plan 4x	4x	S	-	U	S	-
E Plan 10x	10x	S	-	U	S	U
E Plan 40x	40x	S	-	U	S	S
E Plan 60x	60x	S	-	U	-	-
E Plan 100x Oil	100x	S	-	U	-	-
Achromat DL 10x Ph1	10x	S	R	U	U	U
Achromat DL 20x LWD Ph1	20x	S	R	U	U	S
Achromat DL 40x Ph2	40x	S	R	U	U	S
Achromat DL 100x Oil Ph3	100x	S	R	U	-	-

(R: Recommended; S: Suitable; U: Usable; -: Incompatible)

✓ **20x objectives for bright-field microscopy**

Plan 20x objectives can be used for bright-field microscopy with 20x magnification.

### Combinations of objectives and 10x (field number 20) eyepieces

#### Compatibility between optical parts and microscopy methods

Objective	Magnification	Total Magnification (x)	Numerical Aperture	Real Field (mm)	Depth of Focus ( $\mu\text{m}$ )	Resolving Power ( $\mu\text{m}$ )	Working Distance (mm)
E Plan 4x	4x	40	0.1	5	115	2.8	30.00
E Plan 10x	10x	100	0.25	2	18.4	1.1	7.00
E Plan 40x	40x	400	0.65	0.5	2	0.4	0.65
E Plan 60x	60x	600	0.8	0.3	1.2	0.3	0.3
E Plan 100x Oil	100x	1000	1.25	0.2	0.7	0.2	0.23
Achromat DL 10x Ph1	10x	100	0.25	2	18.4	1.1	7.00
Achromat DL 20x LWD Ph1	20x	200	0.40	1	6.1	0.7	3.9
Achromat DL 40x Ph2	40x	400	0.65	0.5	2	0.4	0.65
Achromat DL 100x Oil Ph3	100x	1000	1.25	0.2	0.7	0.2	0.23

## 3.12 Microscope Terminology

### (1) Total Magnification

The total magnification of a microscope is the individual magnifying power of the objective multiplied by that of the eyepiece.

### (2) Numerical Aperture (N.A.)

The numerical aperture is an important factor in determining the efficiency of the condenser and objective. It is represented by the formula:

$$N.A. = n \sin \alpha$$

where  $n$  is the refractive index of the medium (air, immersion oil, etc.) between the objective and the specimen or condenser, and  $\alpha$  is the maximum angle at which light enters or leaves the lens from or to a focused object point on the optical axis.

The larger the numerical aperture the brighter the image and the higher the resolution.

### (3) Resolving power

The ability of an optical system to discriminate between two discrete objects separated by a minute distance. The more minute the distance, the higher the resolving power of the optical system. In relation to the numerical aperture, the resolving power is represented by the following formula:

$$\text{Resolving power} = \frac{\lambda}{2 \times N.A.} \quad \lambda: \text{Wavelength}$$

(The resolving power in the above table is indicated for  $\lambda = 0.55 \mu\text{m}$ .)

### (4) Working Distance (W.D.)

The clearance between the front of the objective and the upper surface of the coverglass, when a specimen image is sharply focused. Generally, the higher the magnifying power of the objective, the shorter the working distance.

### (5) Field Number of the Eyepiece

The diameter in mm of the viewfield observable through the eyepiece. When an eyepiece has an indication of "10x / 20", it means that the magnification is 10x and the field number is 20 for that eyepiece.

### (6) Real Viewfield

The diameter of the region of the specimen that is actually observed with the microscope.  
Real viewfield = field number of eyepiece / magnification of objective

### (7) Depth of Focus

The depth (thickness) of the specimen image in focus, extending above and below the focused image plane. The larger the N.A. of the objective, the shallower the depth of focus.

$$DOF = n \times \left( \frac{\lambda}{2 \times N.A.^2} + \frac{250000 \times \omega}{M \times N.A.} \right)$$

$DOF$  : Depth of focus (object side)

$n$  : Refractive index (according to the medium between the specimen and objective)

If the medium is air: 1

If the medium is immersion oil (Nikon Immersion Oil): 1.518

$\omega$  : Resolving power of the eye (assumed to be  $5' = 0.0014$ )

$M$  : Total magnification

$N.A.$  : Numerical aperture of objective

$\lambda$  : Wavelength (for visual observation, assumed to be  $0.55 \mu\text{m}$ )

### (8) Excitation filter

Of the light illuminating the specimen, an excitation filter allows only the wavelengths of light necessary to emit fluorescence from the specimen to pass through. The term "excitation filter" refers to the excitation of fluorescence in this process.

### (9) Barrier filter

A barrier filter allows only the wavelengths of light consisting of excited fluorescence to pass through. The term "barrier filter" is used because the filter acts as a barrier absorbing non-fluorescent light.

# 4


# Assembly

This chapter explains how to assemble the microscope.

Before assembling or connecting devices, thoroughly read the “Safety Precautions” and “Notes on Handling the Product” at the beginning of this manual, and make sure that you follow all the instructions written therein.

**Tools required for assembly:** Hexagonal wrench (with a 3-mm distance between opposite sides) (One hexagonal wrench is provided with the microscope.)

## 4.1 Configuration Diagrams

## 4.2 Assembly

This chapter explains the microscope assembly procedure for each microscopy method.

Before assembling or connecting devices, thoroughly read the “Safety Precautions” and “Notes on Handling the Product” at the beginning of this manual, and be sure to follow all the instructions written therein.

Tools Required for Assembly: Hexagonal wrench (the distance between opposite sides is 3 mm) (one hexagonal wrench is provided with the microscope)

### 4.2.1 Basic Configuration (for Bright-field Microscopy)


#### (1) Removing the fastening plate

Remove the red fastening plate attached to the top of the stage.

The fastening plate retains movement in the front and back directions of the stage and is attached with two bolts. Remove this fastening plate with the hexagonal wrench provided.

✓ **Keeping the fastening plate**

Be sure to keep the fastening plate with two bolts just in case re-fastening is needed.

#### (2) Changing the direction of the eyepiece tube, removing and attaching the eyepiece tube

When you have purchased the microscope as a set, the eyepiece tube is already installed to the main body.

At the time of observation, turn the binoculars toward you and fix the eyepiece tube. At the time of storage, turn the binoculars backward and fix the eyepiece tube.

When changing the direction of the eyepiece tube, or replacing the eyepiece tube to another one, follow the following procedure.



When installing the eyepiece tube to the main body for the first time, also follow the procedure described in “Installing the eyepiece tube.”

#### Changing the direction of the eyepiece tube

Manually rotate the fixing screw for the eyepiece tube a couple of turns to loosen it while supporting the eyepiece tube by finger. Turn the eyepiece tube in the desired direction and securely tighten the fixing screw for the eyepiece tube while the eyepiece tube fits into the round dovetail mount properly.

✓ **Mark of the direction in which to fix the eyepiece tube for storage**

The microscope has an icon on the arm that indicates the direction in which to fix the eyepiece tube when storing the microscope. This mark indicates that the binoculars of the eyepiece tube is turned backward and fixed.

#### Removing the eyepiece tube

Manually loosen the fixing screw for the eyepiece tube while supporting the eyepiece tube by hand.

Then, loosen and remove the stopper screw for the eyepiece tube by using the supplied hexagonal wrench. Remove the eyepiece tube with care not to drop it.

### Installing the eyepiece tube

Insert the eyepiece tube into the circular dovetail mount while tilting it so that the eyepiece tube fits into it properly. Securely tighten the fixing screw for the eyepiece tube while supporting the eyepiece tube by hand.

Fit the stopper screw for the eyepiece tube, and tighten the screw by using the supplied hexagonal wrench.

✓ **Purpose of the stopper screw for the eyepiece tube**

The stopper screw for the eyepiece tube is for preventing the eyepiece tube from falling when the fixing screw for the eyepiece tube is loosened. The eyepiece tube cannot be fastened with only the stopper screw for the eyepiece tube.

✓ **When screws are attached**


If you purchase the microscope without the eyepiece tube, the fixing screw and stopping screw for the eyepiece tube are fitted into the screw holes and taped to prevent them from falling out. When installing the eyepiece tube, peel off the tape and remove the two attached screws, and then follow the procedure described in "Installing the eyepiece tube."

✓ **When using an eye-level riser**


You can raise the eye point by 25 mm by attaching an eye-level riser between the microscope main body and the eyepiece tube.

While tilting the eye-level riser, fit it into the microscope's circular dovetail mount. Fasten the eye-level riser by using your fingers to tighten the fixing screw for the eyepiece tube. Then, tighten the stopper screw for the eyepiece tube by using the supplied hexagonal wrench.

Attach the eyepiece tube to the round dovetail mount on the upper surface of the eye-level riser. Perform the attachment correctly according to the aforementioned procedure.

Installing the eyepiece tube

Eye-level riser attachment diagram

### (3) Attaching the camera mount

The camera mount (C mount) is attached to the camera port of the trinocular eyepiece tube.

Loosen the knurled screw on the camera port and remove the cap and C mount. Photomicrography becomes available after attaching the C mount to the camera, and then attaching it to the camera port.



The optical path of the trinocular eyepiece tube is fixed at the ratio: binoculars : camera port = 50 : 50.

For details about photomicrography, refer to the instruction manual of the camera.

✓ **Camera mount cap**

Do not lose the cap.

When you do not attach the camera, make sure that the cap is placed on the camera port.

Attach the camera mount


#### (4) Installing the eyepieces

The 10x (Field No. 20) eyepieces are attached to the eyepiece tube before shipment.

When changing to a 15x (Field No. 18) eyepiece (optional), be sure to change both the right and left eyepieces together. The right and left eyepieces should be of the same magnification.

##### ✓ Replacement of eyepieces

To remove the eyepieces, loosen the fixing screws for the eyepieces. The microscope does not come with a hexagonal wrench (the distance between opposite sides is 2.5 mm) to loosen the fixing screws. Contact your nearest Nikon representative if you have any questions.

Installation of the eyepieces

#### (5) Installing the objective


Screw an objective into the revolving nosepiece. Screw the objective firmly to the end.

First, remove the specimen from the stage and lower the stage. Remove the objective while holding it by both hands so that it does not fall, and then install a new objective.

Set the objectives in such an order that the objective magnification increases as the revolving nosepiece is turned clockwise (as viewed from the top of the microscope).

##### ✓ Installation order of the objectives

When you have purchased the microscope as a set, three objectives are already installed to the revolving nosepiece. When adding an objective or replacing the objective with another one, also install it in the magnification order as described above.

Example of order of installing the objectives

#### (6) Installing the condenser

The condenser is attached to the microscope before shipment. The condenser is fastened by a leaf spring, allowing attachment or removal without the use of tools.



When attaching or removing the condenser, raise the stage to its maximum height, and lower the condenser to its minimum height.

The condenser can be removed by pulling it towards you.

To attach the condenser, fit it into the dovetail at the mount, and push the condenser in all the way.

##### ✓ Note on removing the condenser

If you pull out the condenser with the centering screw fully loosened, the inner frame of the condenser mount may be removed together. In this case, put the inner frame back to the dovetail, and then tighten the centering screw. Centering work after reattaching the condenser can be performed without problem.


If the inner frame is accidentally removed

## (7) Connecting the AC adapter and power cord

Plug the power connector of the specified power cord into the AC adapter's inlet for equipment.

For details on the specified power cord, see Chapter 7, "Specifications."


## (8) Other components

### 1) Eyepiece Pointer

The eyepiece pointer (optional) serves as a reference for pointing out the specimen to be observed. Attach the pointer to one of the eyepieces. Rotate and remove the viewfield ring from the eyepiece end.

Attach the eyepiece pointer to the eyepiece and then put back the viewfield ring.


### 2) Filter Receptacle and Filters

Two  $\varnothing 45$ -mm filters can be placed on the upper surface of the field lens. (If the thickness is more than 2.5 mm, only one filter is allowed.)

### 3) Centering telescope

Attach the binoculars when adjusting the annular diaphragm position.

When using the centering telescope, remove one eyepiece from the binoculars and attach the centering telescope in its place.

## 4.2.2 Configuration for Phase Contrast Microscopy

When performing phase contrast microscopy, assemble the microscope in the basic configuration for bright-field microscopy, and then attach the following optical parts to the microscope:

### (1) Phase contrast sliders

There are two kinds of phase contrast sliders: Slider for 10x-40x and slider for 100x). Objectives that can be used are limited to the following:

Supported objectives for phase contrast sliders

Phase contrast slider	Objective	
	PH code	Supported NA
For 10x to 40x	PH1	0.25 to 0.5
	PH2	0.55 to 0.85
For 100x	PH3	0.9 to 1.4

### (2) Phase contrast objective (Ph objective)

Attach an objective with a PH code to the nosepiece. See (5) Installing the objective in “4.2.1 Basic Configuration (for Bright-field Microscopy).”

#### ✓ Combination of optical parts

Confirm in advance that the slider is paired with a supported objective with a supported PH code.

### (3) GIF filter (green interference filter)

When performing phase contrast microscopy, mounting a field lens improves the contrast of the phase contrast image. See (8) Others in “4.2.1 Basic Configuration (for Bright-field Microscopy).”


### (4) Centering telescope

Attach the binoculars when adjusting the annular diaphragm position. See (8) Others in “4.2.1 Basic Configuration (for Bright-field Microscopy).”

## Attaching a Phase Contrast Slider

Insert a slider for phase contrast microscopy into the slot on the side of the condenser.

1. **Set the aperture diaphragm lever of the condenser to the NA 1.25 side.**
2. **Pull the dummy slider out of the slot on the right side of the condenser.**
3. **Remove the knob from the phase contrast slider.**  
The knob is in the form of a screw. Unscrew the knob to remove it.
4. **Insert the slider into the slot on the side of the condenser.**  
Turn the slider so that the centering screw faces the front of the microscope.  
The slider has protrusions to prevent incorrect insertion.
5. **Attach the knob to the slider you attached.**  
Attach the knob to the front of the slider.

Insert the phase contrast slider.

### 4.2.3 Configuration for Simple Polarizing Microscopy

When performing simple polarizing microscopy, assemble the microscope in the basic configuration for bright-field microscopy, and then attach the following optical parts to the microscope:

#### (1) Simple polarization polarizer

Attach the simple polarization polarizer over the field lens.  
Adjust the polarization state by turning the polarizer.


#### (2) Simple polarization analyzer

Attach it to the dovetail at the eyepiece tube mount.  
Remove the eyepiece tube and attach the analyzer to the dovetail on the main body (vibration direction: fixed at the microscope front and back direction).


### Attaching the Simple Polarization Analyzer

Install the simple polarization analyzer on the round dovetail of the eyepiece tube mount.

- 1. Remove the eyepiece tube.**
- 2. Place the analyzer in the cavity in the eyepiece tube mount on the microscope main body.**  
Turn the analyzer so that the mark on the analyzer faces the screw on the rear of the eyepiece tube mount.  
The analyzer faces the microscope front and back direction. (This is also the light's vibration direction.)
- 3. Re-attach the eyepiece tube in its original position and fix it in place.**

Attaching the analyzer


Orientation of the analyzer


### Attaching the Simple Polarization Polarizer

Install the simple polarization polarizer on top of the field lens.

- 1. Place the simple polarization polarizer in the upper part of the field lens unit on the microscope main body.**  
If you turn the polarizer unit so that the mark on the polarizer unit faces the front of the microscope, the polarizer's vibration direction becomes the microscope's right-left direction.
- 2. Tighten the fixing screw to fix the polarizer in place.**  
To fix the polarizer in an open Nicol state, remove the fixing screw, move it to another screw hole, and tighten the screw.

Attaching the polarizer

Mark and fixing screw for the polarizer

## 4.2.4 Configuration for Diascopic Fluorescence Microscopy

When performing diascopic fluorescence microscopy, assemble the microscope in the basic configuration for bright-field microscopy, and then attach the following optical parts to the microscope.

### (1) E2-F-FL Dia-FL Set GFP-B

The E2-F-FL Dia-FL Set GFP-B is composed of the following two parts:

- **EX filter slider**

The EX filter slider is installed in the slot on the side of the condenser. This slider is for switching the empty position for bright-field microscopy and the excitation filter for diascopic fluorescence.

- **BA filter slider**

The BA filter slider is installed in the slot on the front of the nosepiece. When you push this slider in all the way, the barrier filter is brought into the optical path. When you pull this slider out until the click position, the escape position for bright-field microscopy is brought into the optical path.

- **Shielding plate**

Attach the shielding plate to the tip of the arm using two fixing screws. This will prevent strong light including UV light from specimen through the condenser from flashing into your eyes.

### (2) Objective supporting diascopic fluorescence microscopy

The following three types of objectives can be used for diascopic fluorescence microscopy:

- CFI E Plan 4x
- CFI E Plan 10x
- CFI E Plan 40x

- ✓ **Restriction on using objectives**

60x and 100x objectives cannot be used for diascopic fluorescence microscopy.

## Attaching the EX Filter Slider


Attach the EX filter slider (excitation filter) to the microscope. The EX filter is inserted into the slot on the side of the condenser.

1. **Set the aperture diaphragm lever of the condenser to the NA 1.25 side.**
2. **Pull the dummy slider out of the slot on the right side of the condenser.**
3. **Insert the EX filter slider into the slot on the right side of the condenser.**

From the operating knob side, insert the slider so that the text "BF" and the abbreviated name of the microscopy method written on the slider are legible from the front of the microscope.

The slider has protrusions to prevent incorrect insertion, and cannot be inserted from the opposite side. To remove the slider, you must operate the knob.

4. **Confirm that the operating knob is on the left side of the condenser.**


Attaching the EX filter slider

## Removing the EX Filter Slider

To remove the EX filter slider, perform the following procedure:

1. **Set the aperture diaphragm lever of the condenser to the NA 1.25 side.**
2. **While pressing the knob on the EX filter slider, push the slider in from the left side.**
3. **Pull out the EX filter slider from the right side.**


Removing the EX filter slider

## Attaching the BA Filter Slider

Attach the BA filter slider (barrier filter) to the microscope. The BA filter is inserted into the front slot of the nosepiece.

1. **Pull the dummy slider out of the front slot of the nosepiece.**
2. **Insert the BA filter slider into the front slot of the nosepiece, making sure the BA filter slider is facing the right direction.**

During insertion, make sure the tilt of the front of the filter matches the tilt of the nosepiece. The BA filter slider cannot be inserted facing the opposite direction.

Attaching the BA filter slider


## Attaching the Shielding Plate

Attach the shielding plate to the arm end.


1. **Turn the eyepiece tube 90°, and fix it as shown in the figure.**  
(See (2) Changing the direction of the eyepiece tube, removing and attaching the eyepiece tube in “4.2.1 Basic Configuration (for Bright-field Microscopy.)”)
2. **Attach the shielding plate to the arm end as shown in the figure, and tighten two fixing screws using a tool supplied with the filter set.**
3. **Turn the eyepiece tube back to the observation position, and fix it.**

✓ **Fixing screws and tool**

The fixing screws and tool are supplied with the Dia-FL Set GFP-B.

- (1) Turn the eyepiece tube 90°.  
(2) Attach the shielding plate using two screws.

- (3) Turn the eyepiece tube back to the original position.

### Attaching the shielding plate

## 4.2.5 Configuration for Dark-Field Microscopy

When performing dark-field microscopy, assemble the microscope in the basic configuration for bright-field microscopy, and then attach the following optical parts to the microscope:

### (1) Slider for dark-field microscopy

### (2) Objective for dark-field microscopy

Objectives that can be used for dark-field microscopy are limited to the following:

- **For dark-field microscopy:** NA 0.65 or lower
- **For observation using a bright-field position:** NA 1.25 or lower

See (5) Installing the objective in “4.2.1 Basic Configuration (for Bright-field Microscopy).”

#### ✓ **Note on the bright-field position of the dark-field slider**

Objectives that can be used in the bright-field position of the dark-field slider have different conditions than for regular bright-field microscopy.

With objectives of NA 1.0 or higher, vignetting might occur in the periphery, causing the image to darken.

## Attaching the Slider for Dark-Field Microscopy

Insert the slider for dark-field microscopy into the slot on the side of the condenser.

1. **Set the aperture diaphragm lever of the condenser to the NA 1.25 side.**
2. **Pull the dummy slider out of the slot on the right side of the condenser.**
3. **Remove the knob from the slider for dark-field microscopy.**

The knob is in the form of a screw. Unscrew the knob to remove it.


4. **Insert the slider into the slot on the side of the condenser.**

Turn the slider so that the position mark faces the front of the microscope.

The slider has protrusions to prevent incorrect insertion.

5. **Attach the knob to the slider you attached.**

Attach the knob to the front of the slider.

Inserting the slider for dark-field microscopy

## 4.3 Setting Up the Microscope

When you have finished attaching the optical parts to the microscope, you must register information about the objectives mounted on the microscope.

### (1) Turn off the microscope.

Disconnect the AC adapter and turn off the microscope.

### (2) Turn on the microscope while pressing the brightness control knob and illumination on/off button.

The SET UP screen appears.

#### ✓ Information on the objective

The address of the nosepiece that is in the optical path and information on the registered objective such as the objective name and magnification are displayed in the OBJECTIVE field.

Browse and select items on the SET UP screen by using the brightness control knob and illumination on/off button. To move the cursor, turn the brightness control knob. To select an item, press the illumination on/off button.

### (3) Input information about the objectives in each position of the nosepiece.

#### 1. Move the cursor (▶) to "OBJECTIVE" and press the illumination on/off button.

The cursor moves to the displayed name and magnification (for example, [E Plan / Mag: 4x]).

#### 2. Turn the nosepiece to bring the address to be registered into the optical path.

When you turn nosepiece, the registered name and magnification of the new objective is displayed. If you purchased the microscope as a set, the objectives included in the set have already been pre-registered.

#### 3. Turn the brightness control knob to specify the name and magnification of the objective.

Information on standard objectives is registered in the microscope.


If you are using a non-standard objective, select an item where only the magnification is displayed.

For positions where no objective is mounted, specify "None."


#### 4. Repeat steps 2 and 3 to register information about all the desired positions.

#### 5. Press the illumination on/off button.


The cursor moves to "ECO TIME."


Turning on the microscope while pressing the brightness knob and illumination on/off button


SET UP screen

Browse items /values Select

1. Select OBJECTIVE.

2. Turn the nosepiece to select the address.

3. Select the objective using the brightness control knob.

#### Objective setting values

	Objective type	Magnification value		Objective type	Magnification value
1	E Plan	4x	10		1x
2	E Plan	10x	11		2x
3	E Plan	40x	12		4x
4	E Plan	60x	13		10x
5	E Plan Oil	100x	14		20x
6	Achro Ph1	10x	15		40x
7	Achro LWD Ph1	20x	16		50x
8	Achro Ph2	40x	17		60x
9	Achro Oil Oh3	100x	18		100x
			19	None	—

**(4) Set how long ECO mode waits to enter sleep mode.**

1. **Move the cursor (▶) to “ECO TIME” and press the illumination on/off button.**  
The cursor moves to the displayed time period (“[15]min”).
2. **Specify how long to wait (in minutes) until entering sleep mode, and press the illumination on/off button (setting range: 1 to 60).**
3. **Press the illumination on/off button.**  
The cursor moves to “SAVE.”

**(5) To save the setting, select “SAVE” and press the illumination on/off button.**

The setting is saved, and the original screen appears.

If you select “CANCEL”, the changes are discarded and the original screen appears.

```

<SET UP>
CANCEL
OBJECTIVE
[POSITION 1 ]
[E Plan Mag: 4x]
▶ECO TIME [15]min
SAVE

```

Setting time before going to sleep mode

```

<SET UP>
CANCEL
OBJECTIVE
[POSITION 1 ]
[E Plan Mag: 4x]
ECO TIME [15]min
▶SAVE

```

Saving the setting

## 5

## Troubleshooting

Misuse of this product may adversely affect performance, even if this product is properly functional. If any of the following problems occurs, be sure to check the following table for possible causes before requesting service.

If you detect problems that are not listed below or the problem still persists after measures are taken, unplug the power cord of the device and contact your nearest Nikon representative.

## 5.1 Problems With the Microscope Image

### 5.1.1 Common Problems for Each Microscopy

#### Dirt or dust in the field of view

Problems	Check Items	Reference
The specimen is dirty if dirt or dust moves when the specimen is moved on the stage.	<input type="checkbox"/> Clean the specimen.	3.3.1
Objective is dirty if dirt or dust disappears when the objective is switched.	<input type="checkbox"/> Clean the objective.	6.1
The eyepiece is dirty if dirt or dust does not move when the stage or objective is moved.	<input type="checkbox"/> Clean the eyepiece.	6.1
Aperture diaphragm is closed too far.	<input type="checkbox"/> Open properly. (→ Section 3.5)	3.8

#### Poor image quality, low resolution, or contrast too low or too high

Problems	Check Items	Reference
No coverglass is attached to the slide.	<input type="checkbox"/> Attach a coverglass 0.17 mm thick.	—
Coverglass is too thick or too thin.	<input type="checkbox"/> Use a coverglass of the specified thickness (0.17 mm).	—
Slide is upside down.	<input type="checkbox"/> Turn over the slide so that the cover glass faces up.	3.3.1
Lenses or specimen are dirty or dusty.	<input type="checkbox"/> Check dirty or dusty parts, and then clean them as appropriate by reference to the problems described in "Dirt or dust in the field of view."	—
Aperture diaphragm is opened or closed too far.	<input type="checkbox"/> Close or open properly.	3.8
Field diaphragm is opened or closed too far.	<input type="checkbox"/> Close or open properly.	3.9
No immersion oil is applied to the tip of the oil-immersion objective.	<input type="checkbox"/> Apply Nikon immersion oil to the objective.	3.4.2
Nikon immersion oil is not used for oil-immersion observation.	<input type="checkbox"/> Apply Nikon immersion oil to the objective.	3.4.2
Air bubbles in immersion oil.	<input type="checkbox"/> Remove bubbles.	3.4.2
Immersion oil adhering to the tip of the dry type objective.	<input type="checkbox"/> Clean the lens.	6.1
No coverglass is attached to the slide for observation.	<input type="checkbox"/> Attach a coverglass.	—

#### Image too bright

Problems	Check Items	Reference
Setting of the brightness control knob	<input type="checkbox"/> Adjust the brightness by rotating the brightness control knob.	3.2.1

**Image not bright enough**

Problems	Check Items	Reference
Setting of the brightness control knob	<input type="checkbox"/> Adjust the brightness by rotating the brightness control knob.	3.2.1
Aperture diaphragm is closed too far.	<input type="checkbox"/> Open properly.	3.8
Field diaphragm is closed too far.	<input type="checkbox"/> Open properly.	3.9

**Darkness at the periphery, no viewfield seen, or uneven viewfield brightness**

Problems	Check Items	Reference
Revolving nosepiece is not in click-stop position. (the objective is not centered in the optical path).	<input type="checkbox"/> Revolve to click-stop position (swing the objective correctly into the optical path).	3.4.1
Dirt or dust on the lens (condenser, objective, field lens, eyepiece), specimen	<input type="checkbox"/> Clean the lens.	6.1 3.3.1

**Image dark on one side (One side of the viewfield (up, down, right, or left) is not in focus.)**

Problems	Check Items	Reference
Revolving nosepiece is not in click-stop position.	<input type="checkbox"/> Revolve to click-stop position.	3.4.1
Specimen rises from stage surface.	<input type="checkbox"/> Stabilize it using the holder.	—
The thickness of the sample (section) is not uniform and has a tilt.	<input type="checkbox"/> Prepare a sample (section) of uniform thickness.	—

**Image shifts during focus (Becomes asymmetrically defocused when moving the focal point.)**

Problems	Check Items	Reference
Revolving nosepiece is not in click-stop position.	<input type="checkbox"/> Revolve to click-stop position.	3.4.1
Specimen rises from stage surface.	<input type="checkbox"/> Stabilize it using the holder.	—
The thickness of the sample (section) is not uniform and has a tilt.	<input type="checkbox"/> Prepare a sample (section) of uniform thickness.	—

**5.1.2 Problems in Phase Contrast Microscopy****Incomplete field of view or no field of view seen**

Problems	Check Items	Reference
The phase contrast slider is not at the click position (the annular diaphragm is not in the optical path).	<input type="checkbox"/> Move the slider until it clicks.	2.2 4.2.2

**No phase contrast image seen**

Problems	Check Items	Reference
You are not using a phase contrast slider.	<input type="checkbox"/> Mount a phase contrast slider and bring the appropriate position into the optical path.	2.2 4.2.2
You are not using a Ph objective.	<input type="checkbox"/> Bring a Ph objective containing a phase plate into the optical path.	2.2 3.11 4.2.2

**Weak contrast in phase contrast image**

Problems	Check Items	Reference
The Ph annular diaphragm of the phase contrast slider and the Ph code of the objective are incompatible.	<input type="checkbox"/> Bring a Ph annular diaphragm that has the same code as the objective's Ph code into the optical path.	2.2 3.10.1
The Ph annular diaphragm of the phase contrast slider and the objective's phase ring image are not at matching positions.	<input type="checkbox"/> Adjust the position of the annular diaphragm.	2.2 3.10.1
The specimen's phase difference is too big.	<input type="checkbox"/> Make a new specimen with a different mounting agent or specimen thickness.	3.10.1
The specimen's phase difference and the Ph objective type are incompatible.	<input type="checkbox"/> Use an appropriate Ph objective for the specimen.	3.10.1
You are not using a GIF filter.	<input type="checkbox"/> Mount a GIF filter on the field lens.	2.2 3.10.1

**5.1.3 Problems in Simple Polarizing Microscopy****Incomplete field of view or no field of view seen**

Problems	Check Items	Reference
The polarizer has been mounted at a tilt.	<input type="checkbox"/> Install the polarizer so that it lies flat.	2.3 4.2.3

**5.1.4 Problems in Diascopic Fluorescence Microscopy****Incomplete field of view or no field of view seen**

Problems	Check Items	Reference
The EX filter slider is not at the click position.	<input type="checkbox"/> Move the slider until it clicks.	2.4 4.2.4
The BA filter slider is not at the click position.	<input type="checkbox"/> Push in the slider all the way. <input type="checkbox"/> Move the slider until it clicks.	2.4 4.2.4

**No fluorescence image seen**

Problems	Check Items	Reference
The excitation filter is not in the optical path.	<input type="checkbox"/> Push in the EX filter slider from the right side, bringing the GFP position into the optical path.	2.4 4.2.4
The barrier filter is not in the optical path.	<input type="checkbox"/> Push in the BA filter slider all the way, bringing the barrier filter into the optical path.	2.4 4.2.4
Diascopic illumination is not bright enough.	<input type="checkbox"/> Adjust the illumination by turning the brightness control knob close to the maximum position. <input type="checkbox"/> Widen the aperture diaphragm.	3.2.1 3.8

**Fluorescence image is dark (even though the light source is turned on).**

Problems	Check Items	Reference
The diascope illumination is not bright enough.	<input type="checkbox"/> Adjust the illumination by turning the brightness control knob close to the maximum position.	3.2.1
The room is too bright.	<input type="checkbox"/> Dim the ambient light.	3.10.3

**Weak contrast in fluorescence image**

Problems	Check Items	Reference
The objective or cover glass is dirty.	<input type="checkbox"/> Clean the lens.	3.3.1
The immersion oil is emitting fluorescence.	<input type="checkbox"/> Use a low-fluorescence oil (Nikon immersion oil) specified by Nikon.	3.4.2
The slide glass is emitting fluorescence.	<input type="checkbox"/> Use a non-fluorescent glass slide.	—

**5.1.5 Problems in Dark-field Microscopy****Incomplete field of view or no field of view seen**

Problems	Check Items	Reference
The dark-field slider is not at the click position (the annular diaphragm is not in the optical path).	<input type="checkbox"/> Move the slider until it clicks.	2.5 4.2.5
The objective's NA is not in this microscope's supported range for dark-field microscopy.	<input type="checkbox"/> Bring the dark-field objective into the optical path.	3.10.4

**No dark-field image seen**

Problems	Check Items	Reference
You are not using a dark-field slider.	<input type="checkbox"/> Mount a phase contrast slider and bring the appropriate position into the optical path.	2.5 4.2.5
The objective's NA is not in this microscope's supported range for dark-field microscopy.	<input type="checkbox"/> Bring the dark-field objective into the optical path.	3.10.4

## 5.2 Problems in Operation

### Image cannot be focused with high magnification objectives.

Problems	Check Items	Reference
Slide is upside down.	<input type="checkbox"/> Turn over the slide so that the cover glass faces up.	3.3.1
Coverglass is too thick.	<input type="checkbox"/> Use a coverglass of the specified thickness (0.17 mm).	—

### High magnification objective contacts slide when changed over from low magnification.

Problems	Check Items	Reference
Slide is upside down.	<input type="checkbox"/> Turn over the slide so that the cover glass faces up.	3.3.1
Coverglass is too thick.	<input type="checkbox"/> Use a coverglass of the specified thickness (0.17 mm).	—
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust.	3.6.1
Objective is not attached correctly.	<input type="checkbox"/> Screw the objective all the way in.	4.2

### Difference in focal point too large when switching from one objective to another

Problems	Check Items	Reference
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust.	3.6.1
Objective is not attached correctly.	<input type="checkbox"/> Screw the objective all the way in.	4.2

### Image cannot be in focus.

Problems	Check Items	Reference
The stage vertical movement stopper is lower than the focal position.	<input type="checkbox"/> Turn the stage stopper counterclockwise until it reaches the limit.	3.5.3

### Binocular images are not integrated.

Problems	Check Items	Reference
Diopter is not adjusted correctly.	<input type="checkbox"/> Interpupillary distance is not adjusted correctly.	3.6.2
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust.	3.6.1

### Excessive eye fatigue during observation

Problems	Check Items	Reference
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust.	3.6.2
Diopter is not adjusted correctly.	<input type="checkbox"/> Adjust.	3.6.1
Inadequate brightness or illumination.	<input type="checkbox"/> Adjust brightness using the brightness control knob.	3.2.1

## 5.3 Electrical

### Illumination does not turn on when the illumination on/off button is turned on.

Problems	Check Items	Reference
The AC adapter is not connected to the microscope main body.	<input type="checkbox"/> Insert the DC plug of the AC adapter into the DC input connector.	3.1.1
No electrical power.	<input type="checkbox"/> Insert the power cord connected with the AC adapter into an electrical outlet.	3.1.1

### Flickering or unstable illumination brightness

Problems	Check Items	Reference
The AC adapter or power cord is not correctly connected.	<input type="checkbox"/> Connect correctly.	3.1.1

### Brightness of the illumination changes when objectives are switched.

Problems	Check Items	Reference
The light intensity management function is turned on.	<input type="checkbox"/> To prevent the brightness from changing automatically, turn off light intensity management.	3.2.2

### Illumination turns off when left idle (no operations are performed) for a certain period.

Problems	Check Items	Reference
ECO mode is on.	<input type="checkbox"/> Operate controls such as the illumination on/off button, the brightness control knob, and the nosepiece. <input type="checkbox"/> To prevent illumination from turning off automatically, turn off ECO mode. <input type="checkbox"/> In the settings, change how long to wait before entering sleep mode.	3.1.2

### Objective information not shown on the display, or incorrect information displayed

Problems	Check Items	Reference
Information on the objective is not set.	<input type="checkbox"/> Set the information in the SETUP screen.	4.3
Information on the objective is not updated.	<input type="checkbox"/> Set the information in the SETUP screen.	4.3

## 6

## Care and Maintenance

### 6.1 Cleaning the Lenses

- Dust is best removed using a soft brush or gauze.
- More persistent dirt, such as fingerprints, grease and oil, may be removed with clean lens tissue (or soft cotton, gauze) lightly moistened with absolute alcohol (ethyl alcohol or methyl alcohol).
- To clean immersion oil off the oil-immersion type objective, use lens tissue, soft cotton or gauze lightly moistened with petroleum benzine. If petroleum benzine is not available, use methyl alcohol. In this case, you need to wipe 3 to 4 times because the detergency of the methyl alcohol is somewhat weak.
- Cleaning the condenser  
If the condenser was used with oil immersion, remove the condenser from the microscope and clean the lens as necessary.  
Also ensure that the condenser is removed from the microscope and cleaned as necessary during regular cleaning operations.
- Absolute alcohol and petroleum benzine are quite inflammable. Take great care when handling them and when turning the power switch on and off. Be very careful with fire.
- When using absolute alcohol and petroleum benzene, follow the instructions provided by the manufacturer.

### 6.2 Cleaning the Microscope

- We recommend that you use a silicon cloth to clean the microscope.
- For persistent dirt, dampen a piece of gauze with neutral detergent and wipe lightly.
- Do not use organic solvent (alcohol, ether, or thinner). Doing so may result in discoloration of the coating and plastic parts or removal of the printed text.

### 6.3 Disinfecting the Microscope

- We recommend that you use 70% medical alcohol for normal disinfection of the microscope.
- In case of spillage of a sample to the microscope, determine whether the sample is hazardous. If the sample is hazardous, follow the standard procedure of your laboratory.
- Using organic solvent may result in discoloration of the plastic parts.

### 6.4 Storing the Microscope

- When the microscope is not in use, put a vinyl dust cover over the product to protect it from dust, and store it in a dry place where mold is not likely to form.
- Before covering the microscope with the vinyl cover, unplug the power cord.
- We especially recommend that the objectives and eyepieces be kept in a container (such as a desiccator) with desiccant in it.

### 6.5 Periodical Inspections (Charged)

- To maintain the performance of the microscope, periodical inspections and maintenance are recommended.
- For details, contact your nearest Nikon representative.

# 7

## Specifications

### 7.1 Microscopy (Principles)

Use objectives and eyepieces of the microscope to magnify minute cells and tissue optically, and manipulate levers and knobs of the microscope unit to adjust the focus or move the observation point. Then observe or take photographs of the sample fixed on the slide.

#### **Intended use of the product (for medical care)**

This microscope is intended for use in microscopic experiment and diagnostics of cells and tissues at hospitals or by doctors in private practice in the field of pathology, anatomy, and cytology.

The microscopy with diasopic/episcopic illuminations is used to observe a sample fixed on the slide (cells and tissue) as the specimen.

The product is classified as an in-vitro diagnostic medical device. IVD

This product is not intended for use for measurement.

The scale on the focus knob and stage is an indicator to reproduce the position and does not guarantee the value of the thickness or length of a sample measured using this scale.

#### **Intended user**

It is intended for the medical professional and those who work on experimentations in the field of pathology and cytology.

## 7.2 Performance Properties

<b>Model Name</b>	ECLIPSE Si RS Main Body
<b>Optical System</b>	CFI60 (Infinity-corrected CF optical system) Field Number: It depends on the combination. When using Si eyepiece tubes and other optional products supporting a field number of 20: FN 20 When using Ci eyepiece tubes and other optional products supporting a field number of 22: FN 22
<b>Eyepieces</b>	10X, FN 20
<b>Condenser</b>	Supports phase contrast microscopy, contains an aperture diaphragm, and comes with a slot Aperture diaphragm: Ø2 mm to Ø28 mm Centering range: Ø4 mm or more Vertical stroke: 2 mm upward and 13 mm downward Attachment and removal mechanism: One-touch attachment and removal using a leaf spring
<b>Illumination device</b>	Koehler illumination device using a white LED, containing a fly-eye lens Field diaphragm: Ø1.5 mm to Ø18 mm Field lens: Allows two ø45-mm filters (of thickness 2.5 mm or less) to be attached
<b>Focusing Mechanism</b>	Manual single-axis coarse/fine movement handle (double-sided coarse/fine movement) Fine focus knob graduation: 2 µm/graduation Fine focus knob travel: 0.2 mm/rotation Coarse focus knob travel: 37.7 mm/rotation Vertical movable range: 2 mm upward and 13 mm downward (from the focal plane)
<b>Stage</b>	Height of stage upper surface: 135 mm (above the table surface) Stroke: Right and left axis: 76 mm Front and back axis: 52 mm 0.1 mm graduation readable using Vernier
<b>Specimen holder</b>	Supported slide glass thickness: 0.8 mm or more Fastening method: One sided claw grip type Allows up to two 76 mm x 26 mm glass slides to be attached
<b>Revolving Nosepiece</b>	5-hole fixed type, state-detecting type Tilt angle: 22.5° Objective attachment screw: M25
<b>Supported microscopy methods</b>	Bright-field, phase-contrast, simple-polarizing, diascope fluorescence, dark-field

### 7.3 Physical Properties

<b>Model Name</b>	ECLIPSE Si RS Main Body	
<b>External dimensions and mass (reference)</b>	Main body (only): External dimensions: 236 (W) x 284 (H) x 280 (D) mm Mass (reference): Approximately 5.3 kg Set containing binocular eyepiece tube: External dimensions: 236 (W) x 401.5 (H) x 280 (D) mm (excluding protrusion) Mass (reference): Approximately 6.0 kg	
<b>Electrical shock protection class</b>	Class III	
<b>Environmental conditions</b>	Operation:                      Temperature: 0 to +40°C Humidity: 60% RH Max. (at +40°C) (no condensation) Storage and Transport:        Temperature: -20 to +60°C Humidity: 90% RH Max. (no condensation) Altitude:                         2,000m Max. Pollution:                        Degree 2 Overvoltage category:         Category II Indoor use only	

### 7.4 AC adapter

<b>Names</b>	AC adapter	
<b>Manufacturer's name or trade mark, commercial registration number and address</b>	Adapter Technology Co., Ltd. 6F-9 No. 258 Liancheng Rd Zhonghe District New Taipei 235 TAIWAN	
<b>Model identifier</b>	ATS018T-P050	
<b>Input voltage, input current</b>	Universal 100-240 VAC $\pm 10\%$ , single phase, 0.48A Max.	
<b>Input AC frequency</b>	50-60 Hz	
<b>Output voltage</b>	5.0 VDC	
<b>Output current</b>	3.0 A	
<b>Output power</b>	15.0 W	
<b>Average active efficiency</b>	82.6 %	
<b>Efficiency at low load (10%)</b>	79.6 %	
<b>No-load power consumption</b>	>0.04 W	
<b>External dimensions</b>	50 (W) x 33 (H) x 100 (D) mm	
<b>Mass (reference)</b>	Approx. 170 g (without the power cord)	
<b>Electrical shock protection class</b>	Class I	
<b>Environmental conditions</b>	Operation:                      Temperature: 0 to +40°C Humidity: 20 to 80% RH (no condensation) Storage and Transport:        Temperature: -20 to +80°C Humidity: 10 to 90% RH (no condensation) Altitude:                         2,000 m Max. Pollution:                        Degree 2 Overvoltage category:         Category II Indoor use only	
<b>Safety standards</b>	CE marking, GS mark, C-UL-US Listed, PSE mark	

## 7.5 Power Cord

<b>When used in Japan:</b>	PSE approved detachable power cord set 3 conductor grounding Type VCTF 3x0.75mm <sup>2</sup> , 3m long maximum, rated at 125 VAC minimum)
<b>When used in a 100-120 V region, outside Japan:</b>	UL listed detachable power cord set 3 conductor grounding Type SVT, NO.18 AWG, 3m long maximum, rated at 125 VAC minimum)
<b>When used in a 220-240 V region:</b>	EU/EN standard approved detachable power cord set 3 conductor grounding Type H05VV-F 1 mm <sup>2</sup> , 3m long maximum, rated at 250 VAC minimum)


## 7.6 Safety standards

<b>Safety standards</b>	<ul style="list-style-type: none"> <li>• C-UL-US Listed</li> <li>• FCC Part 15 Subpart B Class A            Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.</li> <li>• CAN ICES-003(A) / NMB-003(A)</li> <li>• Australian EMI (AS/NZS CISPR11)</li> </ul> <p>CE Marking</p> <ul style="list-style-type: none"> <li>• IVD Directive            This equipment complies with the emission and immunity requirements of IEC/EN 61326-2-6.</li> <li>• Low Voltage Directive</li> <li>• EMC Directive</li> </ul>
-------------------------	---

## 7.7 External dimensions (unit: mm)

### Binocular eyepiece tube set

### Trinocular eyepiece tube set