# Trial Transfer Catheter

**en**

INSTRUCTIONS FOR USE

**fr**

MODE D'EMPLOI

**de**

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**it**

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**es**

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**pt**

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**nl**

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**sv**

BRUKSANVISNING

**da**

BRUGSANVISNING

**no**

BRUKSANVISNING

**fi**

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**el**

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# Wallace® Trial Transfer Catheter

## INSTRUCTIONS FOR USE

Wallace® trial transfer catheters are closed ended single-use devices provided sterile and should only be used to assess the passage through the cervix prior to embryo transfer.

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

**NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.**

### WARNING:

- THESE PROCEDURES SHOULD ONLY BE PERFORMED BY PERSONS HAVING ADEQUATE TRAINING AND FAMILIARITY WITH THESE TECHNIQUES. CONSULT MEDICAL LITERATURE REGARDING TECHNIQUES, COMPLICATIONS AND HAZARDS PRIOR TO PERFORMANCE OF THESE PROCEDURES.
- TO BE USED BY, OR UNDER THE DIRECTION OF, QUALIFIED PERSONS IN LINE WITH LOCAL GUIDELINES GOVERNING IN VITRO FERTILISATION, IF APPLICABLE.
- THIS CATHETER IS NOT OPEN ENDED AND NOT FOR EMBRYO TRANSFER. ONLY TO BE USED TO TEST PASSAGE OF CATHETER THROUGH THE CERVIX IN A NON-TREATMENT CYCLE.
- DISCARD IF PRODUCT OR PACKAGING IS DAMAGED.

### Description

These instructions are intended for use with sterile Wallace® trial transfer catheters, product codes

| Catheter Type | Length | Code    |
|---------------|--------|---------|
| Classic       | 18cm   | TT1816  |
| Classic       | 23cm   | TT1816N |
| Sure View®*   | 18cm   | CE418   |
| Sure View®*   | 23cm   | CE423   |

\*Sure View® – for enhanced visibility under ultrasound.

The devices consist of closed inner catheter and a detachable outer sheath that is attached to the inner catheter by a Luer compatible adapter. The inner catheter is 16g with a series of 1cm graduations at the proximal end. The inner catheter protrudes from the outer sheath by 5cm and the outer sheath has a series of 1cm graduations at the distal end. The material of the Sure View® inner catheter is modified to contain small bubbles that improve ultrasound image quality.

### Indications for use

To determine whether the cervix is passable for a Wallace® embryo replacement catheter

### Contraindications

The catheter should not be used:

- In the presence of chronic cervical infection
- In the presence of or after recent pelvic inflammatory disease
- For intra-fallopian procedures

### Cautions

Low resolution ultrasound equipment or use of 'Classic' embryo catheters may compromise the visibility of the catheter.

When performing transabdominal ultrasound, ensure that the patient has a full bladder.

### Instructions for Use

Place patient in the lithotomy, dorsal, or left lateral position. The cervix should be exposed with a speculum and gently swabbed with cotton wool moistened with normal saline or medium.

When using transabdominal ultrasound follow normal procedures such that an optimal image of the uterine cavity is achieved.

With the hubs locked in position, advance the catheter so that the inner catheter passes through the external and internal os, into the mid-uterine cavity. It may be necessary to twist the catheter as it negotiates the cervical canal. When using ultrasound manipulate the probe to obtain the best image of the inner catheter.

### WARNING

**The catheter should never be forced against digitally felt resistance while inside the uterine cavity, as forcing the catheter may result in damage to the endometrial tissue and bleeding.**

**The outer sheath should not be advanced further than the internal os, and should certainly never enter the uterine cavity, as this may result in damage to the endometrial tissue and bleeding.**

**The inner catheter should not be pulled back further than the last graduation. Removing it from the outer sheath completely may result in damage to the endometrial tissue and bleeding.**

If slight resistance is experienced, withdraw the catheter and advance the outer sheath until only the tip of the inner catheter is visible. Aligning the most distal marking on the base of the inner catheter with the hub of the outer sheath ensures a smooth radius protrudes. Using thumb and forefinger curve the outer sheath to complement the patient's anatomy and pass the catheter in this position to the internal os. Advance the inner catheter into the mid-uterine cavity. In the rare event where the cervical canal cannot be negotiated by this method, stylets are available to assist insertion. Use code no 1816ST with 18cm catheters and 1816NST with 23cm catheters.

Place tip of catheter approximately 1cm from the fundus. Assess the passage of the catheter in preparation for embryo transfer.

Remove and dispose of the catheter in accordance with local medical hazardous waste practices.

The above are only suggested instructions for use. The final decision on the technique used is the responsibility of the clinician in charge.

\*The products described are covered by one or more of the following: U.S. Patent No. 8,092,390 and counterpart foreign patent(s). Pat pending. Patent applied for.

Sure View, Wallace, Wallace design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark Office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners. © 2014 Smiths Medical. All rights reserved.

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*Sure View, Wallace, designmerket Wallace og designmerket Smiths Medical er varemerker som tilhører Smiths Medical. Symbolet® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive eiere.*