

TPC ICSI Micropipettes

- Micropipettes since 2000
- No compromise on quality and consistency
- Designed for optimal outcome and minimal oocyte degeneration



*"TPC pipettes are some of the best pipettes
I have used."*

Dr. Jorge Cuadros Fernández,
Clinical Embryologist, FIV Madrid, Spain

TPC ICSI Injection Micropipettes

Advantages



- Easy immobilization of sperm
- Very fine fluid control and suction during use
- Easy penetration of the zona pellucida and gentle progression through the oolemma with minimal risk of tearing
- Optimal design for loading and injecting multiple sperm

Characteristics



- Very long taper
- Long parallel wall section
- Long bevel
- Sharp-Blunt Spike™

ICSI Micropipettes

Type	Reference no.	Available angles	Distance tip to bend	Bevel	Inner diameter
Standard Injection	LISR	30	1 mm	30°	5 µm
Standard Injection	LICR-TA35	35	1 mm	30°	5 µm
Standard Injection	LICR-TA25	25	1 mm	30°	5 µm
Standard Injection	LICR-TA20	20	1 mm	30°	5 µm
Small ID Injection	LICR-ID4	30	1 mm	24°	4 µm

The Pipette Company pipettes are CE marked.

Copyright The Pipette Company • Order No. 60010335 • Version 2 : June 13, 2017.

All information correct at time of print. Specifications are subject to change without notice or obligation on the part of the manufacturer.

ORIGIO a/s
Knardrupvej 2
2760 Måløv
Denmark
Tel: +45 46 79 02 00
www.origio.com



CooperSurgical Fertility Companies