



OOCYTE RETRIEVAL BY ORIGIO®

ORIGIO Stiff Single Lumen Needle

Oocyte Recovery System



Kink resistant tubing

Echo markings



Diamond powder bevels
for a sharp finish



Echo markings
at needle tip



Single lumen
handle

CE
0543

origio

a CooperSurgical Fertility Company

Our needles are produced to the highest standards using diamond powder bevels for a sharp, accurate finish. All needles are fitted with kink resistant tubing with needle gauge identification.

- Stiff needle for greater control during procedure
- Improved echo marking for greater visibility under ultrasound
- Triple cut needle tip to prevent excessive bleeding and coring
- Kink resistant tubing with needle gauge identification
- Uniquely designed hub for greater control during oocyte retrieval
- Available in 16g and 17g
- Needle length: 30cm or 35cm

Technical Information	16g	17g
Needle (ID)	1.06mm	0.92mm
Needle (OD)	1.65mm	1.47mm
Volume of needle not including hub	0.32mm ³	0.24mm ³
Volume per cm length of needle	0.009cm ³ /cm	0.007cm ³ /cm
Length of echogenic markings	8mm	8mm
Aspiration tubing (ID)	1.70mm	1.50mm
Aspiration tubing (OD)	2.50mm	2.50mm

Product Compatibility and Order Codes:

	Stiff Wall Type Single Lumen Needle															
	SL1603011S - Quantity: 25 16g Needle 30cm	SL1603511S - Quantity: 25 16g Needle 35cm	SL1603014S - Quantity: 25 16g Needle 30cm	SL1603014LS - Quantity: 25 16g Needle 30cm	SL1603064S - Quantity: 25 16g Needle 30cm	SL1603044S - Quantity: 25 16g Needle 30cm	SL1703011S - Quantity: 25 17g Needle 30cm	SL1703511S - Quantity: 25 17g Needle 35cm	SL1703014S - Quantity: 25 17g Needle 30cm	SL1703064S - Quantity: 25 17g Needle 30cm	SL1703044S - Quantity: 25 17g Needle 30cm	SL1703011MS - Quantity: 20 17g Needle 30cm	SL1703011NS - Quantity: 20 17g Needle 30cm	SL1703014NS - Quantity: 20 17g Needle 30cm	SL1703011BS - Quantity: 20 17g Needle 30cm	SL1703014BS - Quantity: 20 17g Needle 30cm
Fitted with 40cm aspiration tube						✓						✓				
Fitted with 60cm aspiration tube					✓					✓						
Fitted with 100cm aspiration tube	✓	✓	✓	✓			✓	✓	✓			✓	✓	✓	✓	✓
Fitted with 40cm suction tube			✓	✓	✓	✓			✓	✓	✓			✓		✓
Fitted with 100cm suction tube	✓	✓					✓	✓				✓	✓		✓	
Fitted with stopcock in middle												✓				
Fitted with stopcock at needle													✓	✓		
Fitted with stopcock at bung															✓	✓
Fitted with large bung				✓												

✓ = Compatible

Legal Manufacturer: Smiths Medical International Ltd, Boundary Road, Hythe, Kent CT21 6JL, UK

ORIGIO a/s
 Knardrupvej 2
 2760 Måløv
 Denmark
 Tel: +45 46 79 02 00
 www.origio.com



CooperSurgical Fertility Companies